

Product datasheet for SC210217

UBAP2 (NM_018449) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: UBAP2

Synonyms: UBAP-2

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_018449

Insert Size: 814 bp

Insert Sequence: >SC210217 3' UTR clone of NM_018449
The sequence shown below is from the reference sequence of NM_018449. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site **Blue**=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

TACGGCAACTCTCCATACTGGACAACTAAACCCAGAAGAGAGGGGTGGGCTGGGGCAAGGCTTATCCTG
GGCAGGAGAGAACACACGAGCACGTATTTGGGAGCCAGTGCCCTTTCCTAGAAATCCCGACATGTGTCA
GCCATGCCCTCTGTGGGAGTCTGCCTCCAGACTGGCTACTGTATGTAATGTATTTATGTATGATTTGT
AAATGTGATAGAAGTCTGGGGGGAGTTGGGGGATGGCGGCAGATGTTAGCCAGGTCTGCCCTCCCAT
CAAGCCCCCTTCTCACTGTAGCAAAATAAGCACCCCCACCCATCTGCCTTCAGGTCTTCTCACAGCCT
GCACTGCCAGTGGGCCACTAGGGGCAGTCTCTGGAAGGGCTGGTTCAAGGCTGTTTGGGTATAGGGGTC
AGGTACCAATGAAGAATCACGACTTGTCTCACTCCTTTGGAATTGTTTCTTCTGTGTAATTACTTC
ATACCTCTGTTTTGAGAACTGTTCCGTTTGTCTATCTGTCTGCTCCTTCCACCAAAATCTCATCTG
GGAATAGCAGCGGTATCCCTCCACCAAGTATGGCCACCTGTTTGTCTTCATATAGAACAGGGGCTTCTG
GTCTGGTCATGTCCCTAGAGACTTACTAGAGACTGGCTGACCATGAAAAAAAAAATAATCTGAGTTA
TTCAGCCAGTCAGTCCCCTTTAAAAATCTAAGGTATCTCCAAGCAACACATCTCCATTCTAGTTGTCTT
GAACCACGTGGGCAACTGTTTCAGGGGTTTCTCTTGCCCAATC

ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018449.2</u>
Summary:	The protein encoded by this gene contains a UBA (ubiquitin associated) domain, which is characteristic of proteins that function in the ubiquitination pathway. This gene may show increased expression in the adrenal gland and lymphatic tissues. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
Locus ID:	55833