

Product datasheet for SC210226

PET117 (NM_001164811) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PET117 (NM_001164811) Human 3' UTR Clone
Symbol:	PET117
Synonyms:	CSRP2BP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001164811
Insert Size:	842 bp
Insert Sequence:	>SC210226 3'UTR clone of NM_001164811 The sequence shown below is from the reference sequence of NM_001164811. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TTATTGGCAAAAGGATCTCAAAATCATGACTTGAATGTGAAATATCTGTTGGACAGACAACACGAGTT
TGTGTGTGTGTGTTGATGGAGAGTAGCTTAGTAGTATCTTCATCTTTTTTTTGGTCACTGTCCTTTTA
AACTTGATCAAAATAAGGACAGTGGGTCATATAAGTTACTGCTTTCAGGGTCCCTTATATCTGAATAAA
GGAGTGTGGGCGAGACACTTTTTGGAAGAGTCTGTCTGGGTGATCCTGGTAGAAGCCCCATTAGGGTCAC
TGTCAGTGCTTAGGGTTGTTACTGAGAAGCACTGCCGAGCTTGTGAGAAGGAAGGGATGGATAGTAGC
ATCCACCTGAGTAGTCTGATCAGTCGGCATGATGACGAAGCCACGAGAACATCGACCTCAGAAGGACTG
GAGGAAGGTGAAGTGGAGGGAGAGACGCTCCTGATCGTCGAATCCGAGGATCAGGCATCAGTGGACTTA
TCGCACGACCAGAGTGGGGATTCCCTCAACAGTGATGAAGGAGACGTGCTTGGATGGAGGAGCAGCTG
TCCTACTTCTGTGACAAGTGCCAAAAATGGATACCAGCCAGTAAGGAGCTTCTCAATTCCTTTGATTG
TCAATTCCTGTGTGAAGTTTGTTCACCTGTGAAAGAAACGTGAATGTAAAGAGACCTAAATA
AAAGGATAATTATATTTATCTCTAGTTGATCAGCTATAAATTTATATAAAACATAGGCATGTTTGTAC
TAATGAAACGTACTGTCAACCTATCACATTGTAAATTAACACTTTTGGTGGTAACTCAATAAAAT
GAGAAAATTGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001164811.2</u>
Locus ID:	100303755
MW:	31.5