

Product datasheet for **SC210228**

Menin (MEN1) (NM_130803) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Menin (MEN1) (NM_130803) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MEN1
Synonyms: MEAI; SCG2
ACCN: NM_130803
Insert Size: 857 bp
Insert Sequence: >SC210228 3'UTR clone of NM_130803

The sequence shown below is from the reference sequence of NM_130803. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCCTCAAGCGGCAGCGCAAAGGCCTCTGACTACTGGGGACTTCGGACCGCTTGTGGGGACCCAGGCT
CCGCCCTTAGTCCCCCAACTCTGAGCCCATGTTCTGCCCCAGCCCAAAGGGGACAGGCCTCACCTCTA
CCCAAACCTAGGTTCCCGGTCCCGAGTACAGTCTGTATCAAACCCACGATTTTCTCCAGCTCAGAACC
CAGGGCTCTGCCAGTCTGTAATATAGGTCTCTTCTCCAGAAATCCAGCCGGCCAATGGAACCT
CACGCTGGGTCTAATTACAGTCTTTAAAGGCCAGCCCTAGAAACCAAGCTCCTCCTCGAAACCG
CTCACCTAGAGCCAGACCAACGTTACTCAGGGCTCCTCCAGCTTGTAGGAGCTGAGGTTTACCCTTA
ACCAAGGAGCACAGGTCCACCTCCAGCCGGGAGCCTAGGACCACTCAGCCCTAGGAGTATATTTT
CGCACTTCAGAATCCATATCTTGCGAATCCAAGCTCCCTGCCCAAATAACTTCAGTCTGCTCCAGA
ATTTGAAATCCTAGTTTCTCTCTCGTATCCCGAGTCTGGGACACAAAACCTCCGCCCCAGCCTAT
GAGCATCTGAGCCCGCCTCTTCTGACGAACTGGCCCGGATCAGAGCAGGACCTCCCTTCCGAC
CCTCTGGGAACCTCCAGAGTCCAGCCATCTCGGAGCATCCCGGAGGAAATCTGCAGAGGTTAGGA
GTGGGTGACAAGAGCCTGATCTTCTGTTTTGTACATAGATTTATTTTTCAGTTCCAAGAAAGATGA
ATACATTTTGTAAAAAATAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_130803.2</u>
Summary:	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]
Locus ID:	4221
MW:	31.2