

Product datasheet for **SC210232**

Menin (MEN1) (NM_000244) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Menin (MEN1) (NM_000244) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MEN1
Synonyms: MEAI; SCG2
ACCN: NM_000244
Insert Size: 857 bp
Insert Sequence: >SC210232 3'UTR clone of NM_000244

The sequence shown below is from the reference sequence of NM_000244. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCCCTCAAGCGGCAGCGCAAAGGCCTCTGAACTACTGGGGACTTCGGACCGCTTGTGGGGACCCAGGCT
CCGCCCTTAGTCCCCCAACTCTGAGCCCATGTTCTGCCCCAGCCCAAAGGGGACAGGCCTCACCTCTA
CCCAAACCTAGGTTCCCGGTCCCGAGTACAGTCTGTATCAAACCCACGATTTTCTCCAGCTCAGAACC
CAGGGCTCTGCCAGTCGTTAGAATATAGGTCTCTTCTCCAGAAATCCAGCCGGCCAATGGAAACCT
CACGCTGGGTCTAATTACCAGTCTTTAAAGGCCAGCCCTAGAAACCAAGCTCCTCCTCGAAACCG
CTCACCTAGAGCCAGACCAACGTTACTCAGGGCTCCTCCAGCTTGTAGGAGCTGAGGTTTCACCCTTA
ACCAAGGAGCACAGGTCCACCTCCAGCCGGGAGCCTAGGACCACTCAGCCCTAGGAGTATATTTTC
CGCACTTCAGAATCCATATCTTGCGAATCCAAGCTCCCTGCCCAAATAACTTCAGTCTGCTCCAGA
ATTTGAAATCCTAGTTTCTCTCCTTCGTATCCCGAGTCTGGGACACAAAACCTCCGCCCCAGCCTAT
GAGCATCCTGAGCCCGCCTCTTCTGACGAACTGGCCCGGATCAGAGCAGGACCTCCCTTCCGAC
CCTCTGGGAACCTCCAGAGTCCAGCCATCTCGGAGCATCCCGGAGGAAATCTGCAGAGGTTAGGA
GTGGGTGACAAGAGCCTGATCTTCTGTTTTGTACATAGATTTATTTTTCAGTTCCAAGAAAGATGA
ATACATTTTGTAAAAAATAAAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_000244.3
Summary:	This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]
Locus ID:	4221
MW:	31.2