

Product datasheet for **SC210238**

ANKRD34A (NM_001039888) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ANKRD34A

Synonyms: ANKRD34

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_001039888

Insert Size: 858 bp

Insert Sequence: >SC210238 3'UTR clone of NM_001039888
The sequence shown below is from the reference sequence of NM_001039888. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTAGGTGGGCTGGGGAGCCAGGGCGCTGAGAGGAGTGAGAGGAGCCAAGGAGGGTGGGGATTAGGACC
TTAATGGGATAGGTGGGAGAACAGCTCACACACACACCACGAACAAAGGGATCTCTTGCAATGTGCTCA
TGGGTGCACAGTGCACCCACACATGGTTAAACATGTGTCCATAAGCACAGTACTCTTATTAGCAAGAAA
GAATAAGTACTCTATTTACTGGGTATATACACACATGTATTATGCCCTGGTCATTTACATGCATATG
GGATTACTTGTGTACCTACAGGCACAGCATATCAATAAACACACAGGGCTAGGGTCCCAAACTCTCTTG
CATTTGTTCAGAAGCAGTAGGGAACCTTACTTATGGGCAGTCTCCGATGTCTTTGCCCTCCTGCAGAA
GCAATCCCTTGCTTTCCATGTATGCCCTGAAACACAGCCTATCATGATTTGCACACTCTCTCCGTGGA
CATCCAGCCCATCTGCAGGTCTTGCTTCTTTCTCCAGGCCTGGCTCCTTATTCACACCCTGCTTCTA
ACCCCAACCACTTTTCAGGATATCTTTTTCACTCTTCTGGGGTTCCTGCAGCAACCCCAAGCCTCAG
TTCTGCTGCTGCTTTACTGCTCTCAACTTCACTTCTCAACTTCTGCTTTTATCCACCCCTGTTCC
TCCAACAACTACAACAGATTGTCTATTTTCCAGGGGGTGGGTCATGGGCTCTAAGCTGCTCTTCTGT
CTGTCTGAGCTTATGAACACCCCAACAGGTAGAAGTGGGGAGAAACCCTAAATCACTCCTCTCCACT
TGACATATCAAATAAATGTGTTCAAGTC
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001039888.4</u>
Locus ID:	284615
MW:	32.1