

Product datasheet for **SC210260**

TRAFD1 (NM_001143906) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: TRAFD1

Synonyms: FLN29

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_001143906

Insert Size: 843 bp

Insert Sequence: >SC210260 3'UTR clone of NM_001143906
The sequence shown below is from the reference sequence of NM_001143906. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GGGGATGCAGAAGAGGAAGAGGAGGAGTAATGGTGTCTCCAGAGACTTTACATCGGTTCTGTCTTCTG
TGCACAGCAGCACTTGCCGCTGTGCAGGCCACCTCTTTGGCTCTTTGGGTGGGAGAGTTTTCCAGAT
TTTAGATTTTTCTAGGTTATGGCCATTTTGTCTTTTGTGAGTTGTGCTGTGGGGTTTGGGTTTGAGG
GAAGGGAGCAGGGTGGCGGTTGAGGAACGCTTCAGCCTTAGCTGCTACCTTTCGGCAGCAGTGAAATAC
AAGCTGCAGCCTCGGCTGCCAGGGCTCCCTTTTGAATTATTGTCGCCACTGCCCTTGGTGTGTGTGG
TCCAGTGGAAGGAGGGGAAGATTTTGAAACCTGGTAGCCACCAGTAAGGTGATTCTCTGCCCTGTTG
GGGCCATAAATTTGGGGCTTTTGGGCAACCTCTCCGTGTACTGCGTCTGTCCACTCGATTGGGCCCC
AGGTGTGTATGAGGCGCTCTGGTAAGGTGCTCAGGCCAGTTGCAATGTCTGTCAAGAACGAGCTTTTG
ATGTGTTGAGCTGGAGGTGAGTGGACCGGGGCTGTGTTTTAAGCTGCTTCCTTGGCATTGTCATCAC
TGCTTCTGTTCCCGGGGAGCATGGATCTTTTGTCTCACTGCTTTCTAATGGGGAGGGCTGAGGGCT
CCCTGTCCCCACAGCAGGTATGTTGGGCTCTGCCCCAGCCCCCACTTGCTCTGAAAACCAAGTGTGAG
AGCCCCCTCCCTTGTGTTTTTACTGTTATAATAATTATTAACCTTCCTTGAATAGAAATAAAGTT
TGTAATTGGAGTTCA
ACGCGTAAGCGGCCCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001143906.2</u>
Summary:	The innate immune system confers host defense against viral and microbial infection, and TRAFD1 is a negative feedback regulator that controls excessive immune responses (Sanada et al., 2008 [PubMed 18849341]).[supplied by OMIM, Dec 2009]
Locus ID:	10906
MW:	31