

Product datasheet for **SC210269**

CXX1 (FAM127A) (NM_001078171) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CXX1
Synonyms:	CXX1; FAM127A; Mar8; MAR8C; Mart8; MART8C; SIRH5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001078171
Insert Size:	851 bp
Insert Sequence:	>SC210269 3'UTR clone of NM_001078171 The sequence shown below is from the reference sequence of NM_001078171. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTTGGATGGGAGGAGGACGAGGACTTCTAGGCCGGGAGACCCCTCGGGCCTGGGGGCGGGTGTCTGGGG AGGGTCCGCTGTGTTACTGGCCGCCGCCAGGGTCGCCACCGGCCCTCCCTCCGCGCTCCCTCCCCC TCGAGCCGCCGCGATGTCCCTGCGCTCCTGTTCCCTCCCGCTAGTGCTTGCTTTGTTCCAGGAATA GCGCTCCAGGCTCCTGCTGCCGCCCTGGGCTCACTCTGGAGCGAGCCGCCCTCTCCTTCCAGCC AGCCAGCCCCCTCCATGTACATTTGGACGCTGTCTGCGCTCCAGCTGCAAGCTGGGCTCCTGTTACAC ACTGGACAGACCACCACTGCCGCCCTGCCAAGCCCTCTCTCCCAACAGACTGCCAGACGACTACA TCATTCTGCCACAGACTGCGCTGCCACAGCCATCGCCATCCATCGCATCCCAACGACAGACTGCTGC TCCTAGTGATCTGGAATCACCTCGGAGGTATCTGGGCTGGCCACAGTCCCTGGACAGTGATCCAGACAG CTGGCCGCCCCCCAAGGGATCTGTACCTTCAGCGAGACCTATTTCTCCCAACCCCAAGAACTCTT GTGTTCTTGCTAGGCCAGGTGTTCTGGCAGCCAAATCGAGTCTCTCATTTTCTCTTGTGGACAGT TAGTTTTGCCATAACGAGTATTCTGAGTTTGCAACTGTCTCTGATGTGTGCTTTTGTCAACAC AGTAACCCCTGCATTCTGCTCTGCTCTAATACTACCTACCTGGAGAAAGTCTTTCTTATTTCAATAAA TGTCAGACATTATTGAAAAGAAA ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001078171.2</u>
Locus ID:	8933
MW:	30.5