

Product datasheet for SC210287

ACAD11 (NM_032169) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ACAD11 (NM_032169) Human 3' UTR Clone
Symbol:	ACAD11
Synonyms:	ACAD-11
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032169
Insert Size:	860 bp
Insert Sequence:	>SC210287 3'UTR clone of NM_032169 The sequence shown below is from the reference sequence of NM_032169. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAAGCCAAAAGACTGACAGCCAAGATATAAGGAGGGTGGCACTGCCACATCCCACTGGCAGAACTCTC
CTTTATACAACTTCATTGGCTCCAACATTTGAATCTCATATTTTGTAGCAGTTTGAGCACAGGGTTA
ATTATTCATTTGTGGTAAAGATTATAGCATCTATTTTGATCAGTGGGTTTTATTATTTCAAGGGTCACA
CAGGGTTAAGTTCAGTAAGAAATGCTGTAGCTGTTGTCATTCAATCTAGTGCCTCCTTGAGGCCAGGAG
TTCAGGACCAGCCTGGGCAACATAGCGAGACCCCATCTGCTACAAAAATTTAAAAATGAAACAAGTGT
GGTGGCACATGCTTGTAGTCCTAGCTACTTGGGAGGCTGAGGCAGGAGGATTGCTTGAGTCTAGGAGTT
TGAGGTTACAGTAAGCTGTGATCGTGACATGGCCTCCAGCCTGGGTGACCGAGTGAGACTGTTTCTAAA
AATAAAAAACAAAAATAAATTTCTTCTTGAGGTGGGTGGAGGTGGGGAGCAAGAATTTGACCTGGCTC
TGATCCCTGGTGTGTGTGTGGCCTCTTAACGTTTGCCACTGAGCCTTAACCTCACTGTACTTCACT
GTACTTCACACGATTGGTGTTAACATTTTAATCTTAGAAGACCTGACCACTGAGGGTTTGTGTGA
GAATTGCTGAAGCCACGTAGAAGCACCTTGAAATCTGTAAACCACAAGAAAGTACTTTATAAAAGGTA
TCCTTATTTGAAGTGGATAAATCTGTAACTCGAAAAGTTGTGATTTAGAAGACAGGATTGTTTTGAA
CATTAGGAATTAAGGCTATATCTGGTCCTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032169.5
Summary:	This gene encodes an acyl-CoA dehydrogenase enzyme with a preference for carbon chain lengths between 20 and 26. Naturally occurring read-through transcription occurs between the upstream gene NPHP3 (nephronophthisis 3 (adolescent)) and this gene. [provided by RefSeq, Aug 2015]
Locus ID:	84129
MW:	31.9