

Product datasheet for **SC210309**

TAGAP (NM_152133) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TAGAP (NM_152133) Human 3' UTR Clone
Symbol:	TAGAP
Synonyms:	ARHGAP47; FKSG15; IDDM21; TAGAP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152133
Insert Size:	1376 bp



[View online »](#)

Insert Sequence: >SC210309 3'UTR clone of NM_152133
The sequence shown below is from the reference sequence of NM_152133. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCCAATATGCCAAAGAATCGTATATTAGGAGGGAGGCCATACGCCATGCCATAGCTTGTGCTATCTG
TAAATATGAGACTTGTAAGAAGTGCCTGTAGATTGTTTTTAAAAGGTCTTGAATAAGCTCCTTGAGAA
AGTTGTGGAAAGCCCTCCTCAGTGAGGATAGCTACACCATGGCCATGGCGCATCAGATAGTCTCTGTGT
ACCTGGATTTGTGCAATATGTAAAAATGTATCAAAATGATTATAGATAAGGTGTTAGGTGCAAAGGATG
TCTAATAATCCCTGCACACGTTTTGAAGTGCAGTGAAGTACACTGCTGTTCTTCTGCTTCTGGGCAC
TTTTCTCTTGGTTAGTGTAAAAATTATCTTCGCTTTTTTAATGTGGCCTCAAATGTCATGCCAATTT
TCACATCTTCCACAACTCCATTTAGGGAGAAATGTTAAATCTCTGGTATAAGTTTACTCCATACCAG
AGTAAACTATATATTACTCTATATAAGCAGTCTTGCAATAACTAATCACCACCATAGAAGAAAGAAACA
GACTGCAAGGAACAGAGTTGAGTGTCTGGAGTCATCAAAGGCATTAATAAACTCCAGTAAAGCTGGGGC
CGTAGCAAAAATCATGAAAAACACTTCAACGTGTCCTTTCAATCATCCAATTAATGTGGGTAGATTAA
TGAAAAATGTATTACATCAATATTAACCTCATCTATAGCACTTTGAGTATCTTTGTAGTTTCATGATCCT
ATCCTATAATGTGGAGGTAAATGATTTTATATGCATTGGGGGTCAATATAAACTTCAATGTAATTTTC
ACTACAATAAATTGCCTTCCTTATTTGAAAGTAAAAATGTTCTTTTTTTATTGACATGTTAACTCCTTT
CCCTCAATTAATAGAAATCCACTAAGAAATCATTACATATTGGTTCACTCAACAAGCATTATATAAA
TATATATTCACATATTCTAGACTAATAGCAAGACCGGGGATCTTGTTTAGGAAGAAGCAGTCCCTGTTCT
CCAGAATTTGCAAAACCATTAATAAAGCACCTACTTTAAGCCATTTTTTTCAGCAAGGAGTCATTCTG
CCAGAAAAATGTAGTACACAAATACAGGATAATAACAAATGTAAAAATTCTCATTTCTAGTGAATTA
AACTTTCCAGTAATTTTACATCTCACCTCATTCTTATGATGCTCAGTTTGCTTAATTATTGGCAACAT
AATGGTAAAAATGTTGTACTGTATTAGAGCTTTACTGTTGATTATTAAGATATTTATCCAGTATCTTA
GAGATGCTGGACTTCAATTTTCTTATTTTATCTTTTTAAATAAAATATCATGCATGCTTCTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_152133.3](#)

Summary: This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013]

Locus ID: 117289

MW: 53.3