

Product datasheet for **SC210316**

ZNF540 (NM_152606) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF540

Synonyms: Nbla10512

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_152606

Insert Size: 946 bp

Insert Sequence: >SC210316 3'UTR clone of NM_152606
The sequence shown below is from the reference sequence of NM_152606. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CAACATCAGAAAACTCATAATGTAATTAAATATAAGAAAAGGTTTCCATGTCATGCTCTATTTATAGAA
TATCAAAATATTTATGGCCAGAAGTTCTGTCAATGTGTTGATGTTTTTACACATATTAACCTAATAA
ATGTATGAGTCTTAAATACCTCTTAGTTCTCATTAAATTTAGGAAAATTCACACTAGAAAAATAAGCTG
TTAATGTAACAGTTGTGGAAAAGTGTCTAGCAACAGCATATACTTATCATCATTGCCTTTCCACTACT
CTACTATCTGTGTGATATTAGACAAAATATTTGCTTCTTGGTACCTCAGCTGAAAATGAAACACACCT
AAAAGTGTGGTTGTTTCCAACATGTATAATACAGCAACAATCTGGCCCAAAGTCTTTGGATTAAT
ATTGGATATTACTGTTTTTATTATCATCAACATTATTATTAGTGGATTTCTTAATAGGAAGATGCAATG
GAGATGACAAATTTGGAAAAACCACTCATCACTTACATTTTCAAGTACTTCTTTGATAAAATCTGTT
ATGGGCTGAATGTTGTGTTCCCGTAACAATTCCTATGTTGAAACACGAATCCCAAGGTGATGGTATTT
GAAGGTAGGGCCTTTAGGAGGAAATAGGTCATGAGGGTGGAGCCTTCATGAGTGGAATTACTGCCTTT
ATAAGAAGAAGCCAAAGAGCCAGCTAGCTCTTTCAACCACATGAGGTTACAGCAAGAAGTCAGCAGTCT
ACAGTGCAAAAGAGGGCCTTACCAGAACCAAGCATGCTGGCACCTTGACCTTGGAATTTCAAGCTCC
AAAAGTGTGAGAAATAAACCTCAGTTGTTTATAAGCCACTGAGTCTACGATATTTGTTATGCAGCCCA
GATAACTAAGGCACTCCATTAGAATATAAATAAAATGTAAGATCTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_152606.5</u>
Summary:	May act as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]
Locus ID:	163255
MW:	36.5