

Product datasheet for SC210334

CCK4 (PTK7) (NM_152881) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCK4 (PTK7) (NM_152881) Human 3' UTR Clone
Symbol:	CCK4
Synonyms:	CCK-4; CCK4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152881
Insert Size:	864 bp
Insert Sequence:	>SC210334 3'UTR clone of NM_152881 The sequence shown below is from the reference sequence of NM_152881. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
GGAGACAGCACCGTGGACAGCAAGCCGTGAGGAGGGAGCCCGCTCAGGATGGCCTGGGCAGGGGAGGAC
ATCTCTAGAGGGAAGCTCACAGCATGATGGCAAGATCCCTGTCCTCCTGGGCCCTGAGGCCCTGCC
TAGTGCAACAGGCATTGCTGAGGTCTGAGCAGGGCCTGGCCTTTCCTCCTCCTCACCTCATCCTT
TGGGAGGCTGACTTGGACCCAACTGGGCGACTAGGGCTTTGAGCTGGGCAGTTTTCCCTGCCACCTCT
TCCTCTATCAGGGACAGTGTGGGTGCCACAGGTAACCCCAATTTCTGGCCTTCAACTTCTCCCTTGAC
CGGGTCCAACCTGCACTCATCTGCCAACTTTGCCTGGGGAGGGCTAGGCTTGGGATGAGCTGGGTTT
GTGGGGAGTTCCTTAATATTCTCAAGTTCTGGGCACACAGGGTTAATGAGTCTCTTGGCCCACTGGTCC
CACTTGGGGTCTAGACCAGGATTATAGAGGACACAGCAAGTGAAGTCTCCCACTCTGGGCTTGTGCA
CACTGACCCAGACCCACGTCTTCCCACTTCTCCTTTCCTCATCCTAAGTGCCTGGCAGATGAAG
GAGTTTTCAGGAGCTTTTGACACTATATAAACCGCCCTTTTGTATGCACCACGGCGGCTTTTATATG
TAATTGCAGCGTGGGGTGGGTGGCATGGGAGTGGGGTGGCCCTGGAGATGAGGAGGGTGGCCAT
CCTTACCCACACTTTTATTGTTGTCGTTTTTTGTTTTGTTTTGTTTTGTTTTGTTTTGTTTTGTTTTA
CACTCGCTGCTCTCAATAAATAAGCCTTTTTTACAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152881.4</u>
Summary:	This gene encodes a member of the receptor protein tyrosine kinase family of proteins that transduce extracellular signals across the cell membrane. The encoded protein lacks detectable catalytic tyrosine kinase activity, is involved in the Wnt signaling pathway and plays a role in multiple cellular processes including polarity and adhesion. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]
Locus ID:	5754
MW:	31.6