

Product datasheet for **SC210341**

TIP120B (CAND2) (NM_001162499) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TIP120B
Synonyms:	TIP120B; Tp120b
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001162499
Insert Size:	851 bp
Insert Sequence:	>SC210341 3'UTR clone of NM_001162499 The sequence shown below is from the reference sequence of NM_001162499. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CCCAGCACAGACTCAATGGAGCTCAGCTAGTCCCCTCAGCACCAAGGTGGGCCCTCGCTTAAGAGAAAG GAGCCCAACCAAGTCCGAGGCCTCCCCATCCCACCATCGCAGGTCTCTACTTTTGCCCTTCACCATCT CACTGGGGGCCCTGTGCTCCTGGTCAGGGCTTACAGTGCCTTCTCCAGGGACCAACTCAAAGGCCCC CAGCCCAAGCTGTGAGGCTGCCAAGTTGGGCCCTTCCCTTAAGTACAGGACAGTCATCCAAAGAAATA GGGTGAGGAAGTTTCCAGTGACTTCACACTGTACCCCTCCATAGTCTGTCTGGTTCCTTCAGAGGGTG TCTCTGCCTCACAAGTAGTAGTATTTAGAAATAGGCTGTGCTGTCAGCTGTAAAGATCAGGAGGCAG CAGACCACTCTGGTTTCTTCACTGCATTAGCAATGCCTGAAGTTAGTGCTCAGGCCGGGCATCTCA AAAGAAAAGATACTTGAGTTATTCATTTTAAATTTAAACGGTTTCATTTTAAAGTGGCAGTGATGA ATCAGAAATTTGGAAGATGATACGGGTTTCTTTTCCAGGGAGGAGGAATGGGTTGGGTAGGGAAGT GACAGGCTTGACCTCATGTTTCAATTTCTAATTTCAAATACTTATTAGCAAATTTGGGCAACAATGGGC ATCTTCCATGCCACCCAGGCATAACCAAGTTGGTTTGTTCCTTCTGAGGAAGGTTTCAATGTGTC TAGTGTTCAGTATTGAGGACAAAGAAATACAAGTGGCAGGCCCAAGTATTTCTGTGATATCCAGGTT AATAAGATTAGATTCTAAGTTA ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001162499.2</u>
Summary:	Probable assembly factor of SCF (SKP1-CUL1-F-box protein) E3 ubiquitin ligase complexes that promotes the exchange of the substrate-recognition F-box subunit in SCF complexes, thereby playing a key role in the cellular repertoire of SCF complexes.[UniProtKB/Swiss-Prot Function]
Locus ID:	23066
MW:	31.9