

Product datasheet for **SC210357**

KCTD1 (NM_001136205) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: KCTD1 (NM_001136205) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: KCTD1
Synonyms: C18orf5
ACCN: NM_001136205
Insert Size: 870 bp
Insert Sequence: >SC210357 3'UTR clone of NM_001136205
The sequence shown below is from the reference sequence of NM_001136205. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCCGGATAAAGCAAGAGCCTCTGGACTAATGGACATATTTCTTATGCAAAAAGGAAAACACACACAA
CCAATAACTCAAACAAAAAGGGACATTTATGTGCAGTTGGGACAGCAAACCAAGTCCTGGACGTAAAA
TCGAATAAAAGACACATTTATATCCAATAGAGACCACACCTGTATTCATATGGGAACAATTGGAATAGT
GATATCCTCAAGGTGTAATAAATATATAAATATATATATATATGTCAAAAGGTAGGAAATGCAAAAAAG
AAAAAAAAAAAAAGGTGACAGCCGAGTTGGTGCTGTGATGGCCGTGAAGTGTCTGGGCCTCCCGAGGC
CTCTGACAAATAAACAAGCCATGAGTGGTGAAGACACAGTCTCCTTACAGTTTCCATTGCCAACACAG
CCATCCATATTTCTTTTTCTTTTGTCTTTCTTTTTCTTTTTTTTTAAAAAACAAAACAAAACAAAAC
ACCTTGAATCAAGTTTGTGTATATGGAGTTCCACGTCTTTCTTTAGGCAGGGACCAGGCAGGACTT
CAGAAAAACCTCATGAGCACATTGCAAAGATGTTAGACATGAAATTTAAATGTAGTTTGTACAGAAG
TCACACTTTTTGTCCACCTCACAGATGTGAACTTACTTTGTTTTAAACTGATCAGTTTTGCCAAGG
GGCCAGAATTATTCCTTGTTAGAAATGCTCCAGTTCAGTCTGCTTTCTTACAATTTTTCAAATTT
TATAATGTATTAATAACAATAAACTCTGTTTTAAAAAATAAGGTCTGTGTGAAACACACATGTGGGGTG
AGGCTGGATTAAGTGAATTTTTCTTTTTGAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001136205.2
Summary:	This gene encodes a protein containing a BTB (Broad-complex, tramtrack and bric a brac), also known as a POZ (POxvirus and zinc finger) protein-protein interaction domain. The encoded protein negatively regulates the AP-2 family of transcription factors and the Wnt signaling pathway. A mechanism for the modulation of Wnt signaling has been proposed in which the encoded protein enhances ubiquitination and degradation of the beta-catenin protein. Mutations in this gene have been identified in Scalp-ear-nipple (SEN) syndrome. [provided by RefSeq, May 2017]
Locus ID:	284252
MW:	34