

Product datasheet for **SC210361**

RNF19B (NM_001127361) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	RNF19B
Synonyms:	IBRDC3; NKLAM
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001127361
Insert Size:	861 bp
Insert Sequence:	>SC210361 3' UTR clone of NM_001127361 The sequence shown below is from the reference sequence of NM_001127361. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

CAGTTCCTACAACCCACAGGACAGGTTTAGCATGATCCATGCA**TGAC**TCAGCAAAGTGGATTTTGTCTCC
ACAGAGAATGCAACAATATGGAATCCAAGTGGACATTGAAGCCAAACCAAGCCACTATCAGCTGGTGAG
TGGAAGCAGCACGGAGGACTCGCTCCATGTTTCATGCTCAGATGGCAGAGAATGAAGAAGAAGGTAGTGGT
GGCGGAGGCAGTGAAGAGGATCCCCCTGCAGACACCAAGCTGTGAACAGAAAGACTGCCTGGCCAGCA
AACCTTGGGACATCAGCCTGGCCAGCCTGAAAGCATCCGCAGTGACCTAGAGAGTTCTGATGCACAGTC
AGACGATGTGCCAGACATCACCTCAGATGAGTGTGGCTCCCCCGCTCCCATCTGCAGCCTGCCCTCG
ACCCCGAGAGCCCAAGGTGCACCGAGCCCAAGTGCCCATATGAACCTCTCTGCCCTAGCCGAGGGACAAA
CTGTCTTGAAGCCAGAAGGTGGAGAAGCCAGAGTATGAAGTGAATGAATGCTCCTGTTCTGAGAAGCAC
ACTTGTAACATGCATCTTTTGAATTTTTTTTTTTTTTTTTTCCAAGGGGTAGAGATTTATGATTTTATTT
CACAGATTCTCTGGTCACAGGTTTTTGCCAGGGAAATCTGAGAAATTCACAAATTTCTTACCAGATAAA
ACATGAAAAGTTTGCCGTTAGTTCCTCCCTCCCTCCCTCCCTCTTTTATGTTTAAATTTATTGGTTAAAC
TGATGGCAGCAATCCATGAGGTGTGTCAAAGAGTGTACATATGTATGTGTATATTGAATGCTAAACAT
ATTACTGAAAGACACATTTTA

ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001127361.1</u>
Summary:	This gene encodes a multi-pass membrane protein containing two RING-type and one IBR-type zinc finger motifs. The encoded protin is an E3 ubiquitin-protein ligase that plays a role in the cytotoxic effects of natural killer (NK) cells. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes X and Y in a possible pseudoautosomal region. [provided by RefSeq, Jul 2014]
Locus ID:	127544