

Product datasheet for **SC210380**

BBS9 (NM_198428) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: BBS9

Synonyms: B1; C18; D1; PTHB1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_198428

Insert Size: 873 bp

Insert Sequence: >SC210380 3'UTR clone of NM_198428
The sequence shown below is from the reference sequence of NM_198428. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GTTTCACCCCTCCAAGGAGTCTCGGAAATTCAAGTAGAGTTGTTGGTTGAGAGGAACATCCCCATC
TCAAGGCCGAACCTGTGTGAACCTCATGCCAAGCACAGATATAGGGCTGGCGCAGGTGCTTCCTAAAGC
TCACCTTCCTGGAGATGACATGCATAGAAAGAGGGTTGGGACTTTTACTTCACTAGGAGAAGTTGTA
ACACCATGGGGAAGTCAGCTGAACTTGTCTGTTTTGCCAGGAAAGGAAGTAGTTGCCTTTGGTCATC
CATCTGCTAATAGTCACAGAATACAGTGAAATGACATAGTTTTGGGTTAGATTTTATAATGCAAAGATT
CAGATCCAAAATAATTTACATCCCATTTTTTACAGAACTTCTATATAGTAAATGTATCAAGTTTAAT
AAAGCATCTCATTGTCAAATAATATCTGGATTTTATTATAATTAGAGGGATTATGAGTGATTGCTC
TACATTATTTCTTCAAAGGAAAGGAAAGGAATTGAAGACTTTGCTACTCTCTGGTAAGACTTGAATGTG
ATTATTTTATAAAATAAAGAACCACATGAAACTTTCATGTGATTGTCTATTAAGCTAAACGTGGAA
GAAGGGAAGTTTTCTAGTATTGTATGAAATGTTAAATATTTTCAAGGAAGCTAGTACATTAACATCC
CCTTTTCATTACCCCATTTGGGTGGCACCTTTTAAAAATCCAGGTTATTAAGGATTACACTTTATTAC
AGTTAAATTTAAGGTTTTTATGTGTAAGAAGTGAAGAAACTTTGCTAATTTCAAGTTCGAATTACAG
TTGTGAGATTGATTTATTTAAATAAAAAATATGCCTTTGCCAGTTT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_198428.3</u>
Summary:	This gene is downregulated by parathyroid hormone in osteoblastic cells, and therefore is thought to be involved in parathyroid hormone action in bones. The exact function of this gene has not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2017]
Locus ID:	27241
MW:	33.9