

Product datasheet for **SC210381**

BBS9 (NM_001033605) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	BBS9
Synonyms:	B1; C18; D1; PTHB1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001033605
Insert Size:	873 bp
Insert Sequence:	>SC210381 3'UTR clone of NM_001033605 The sequence shown below is from the reference sequence of NM_001033605. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GTTTCACCCCTCCAAGGAGTCTCGGAATTAATTCAGTAGAGTTGTTGGTTGAGAGGAACATCCCCATC TCAAGGCCGAACCTGTGTGAACCTCATGCCAAGCACAGATATAGGGCTGGCGCAGGTGCTTCCTAAAGC TCACCTTCCTGGAGATGACATGCATAGAAAGAGGGTTGGGACTTTTACTTCACTAGGAGAATTGTA ACACCATGGGGAAGTCAGCTGAACTTGTCTTGTGTTTCCAGGAAAGGAAGTAGTTGCCTTTGGTCATC CATCTGCTAATAGTCACAGAATACAGTGAAATGACATAGTTTTGGGTTAGATTTTATAATGCAAAGATT CAGATCCAAAATAATTTACATCCCATTTTTTACAGAACTTATATAGTAAATGTATCAAGTTTAAT AAAGCATCTCATTGTCAAATAATATCTTGATTTTATTATAATTAGAGGGATTATGAGTGATTGCTC TACATTATTTCTTCAAAGGAAAGGAAAGGAATTGAAGACTTTGCTACTCTCTGGTAAGACTTGAATGTG ATTATTTTATAAAATAAAGAACCACTATGAACTTTCATGTGATTTGTTCTATTAAGCTAAACGTGGAA GAAGGGAAGTTTCTAGTATTGTATGAAATGTTAAATATTTCAAGGAAGCTAGTACATTAACATCC CCTTTTCATTACCCCATTTGGGTGGCACCTTTTAAAAATCCAGGTTATTAAGGATTACACTTTATTAC AGTTAAATTTAAGGTTTTTATGTGTAAGAAGTGAAGAAACTTTGCTAATTTCAAGTTCGAATTACAG TTGTGAGATTGATTTATTTAAATAAAATATGCCTTTGCCAGTTT ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001033605.2</u>
Summary:	This gene is downregulated by parathyroid hormone in osteoblastic cells, and therefore is thought to be involved in parathyroid hormone action in bones. The exact function of this gene has not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2017]
Locus ID:	27241
MW:	33.9