

Product datasheet for **SC210392**

SLC35E4 (NM_001001479) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	SLC35E4
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001479
Insert Size:	892 bp
Insert Sequence:	>SC210392 3'UTR clone of NM_001001479 The sequence shown below is from the reference sequence of NM_001001479. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CGGAGGGACCAGCCAGCAAGGGTCTTGAGACCTGGGGGATCTCAGGAGCCACCTGGGATGGCCCTGG
CCTGAATCCAGCCTCCGCTGTGGCCATAGAAGGAATGGAGAACAGGGCTGGGCATGGTGGCTCACGCCT
ATAATCCCAGCACTTCCAGAGTCCGAGGTGGTGGATCACCTGAGGCCAGGAGTTCGAGACCAGCCTGG
CTAACATGGCAAAACCTCATCTCTACTAAAAATAGAAAAATTAGCTGGGCATGGTGGCGCTGCCTATA
GTCCCAGCTACATGGGAGGCTAAGGTGGGAGGATCACTTGAGCCCTGGAGATCGAGGCTGCAGTAAGCC
AAGATCGCATGCTACTGCACTCCAGCCTGGGAGACAGAGCGAGACGCTGTCTCAATTAACAAAAAAAAA
AAGTGGAGAACTGGCAGTGACCTCTACTGGGGGCCATGGCAGGGAGGGAGCCTTCTGGAAGGGCTGCC
TTGGAGATTGGAATGGGACTCCCAGGGAGACCTGCGTTCCATCCCTGCCTGCCTCACCCCTGCCACAG
ACTCTGCACACCACTGGATGGTGGGTCCAAGCCTGGCACAGTCCCTGTGCTTGTGAGAGTCATTATTAT
GATTAATATCAATTACGATGCCAAAAATTGCTGGGCAAACTTTGAAGACCTCAACTTGTACAATGACG
ATGATGATGATTCTTGGCGTTACACAATCCTTCTCTGGGGGGGAGGCAGCTAGGAGGCCAGCAGG
GGGGCTTCTATGCTGCTGGGCTCCCTAGGGAGTTGGGGTAGTCTGTGCCAACTCCAGGCAGCTGCTGT
GGCCTCACCCCTGGGCCCCCAATTTTGGGTCATCCATCCTCAAATACACTATTTTTGCTTGTA
ACGCGTAAGCGGCCCGGCATCTAGATTTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001001479.4</u>
Summary:	Putative transporter.[UniProtKB/Swiss-Prot Function]
Locus ID:	339665
MW:	33.1