

Product datasheet for **SC210397**

GRHL3 (NM_021180) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: GRHL3 (NM_021180) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: GRHL3
Synonyms: SOM; TFCP2L4; VWS2
ACCN: NM_021180
Insert Size: 860 bp
Insert Sequence: >SC210397 3'UTR clone of NM_021180
The sequence shown below is from the reference sequence of NM_021180. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAATTCAGATCATCCTTAAGGAGCTGTAGGCCTCTCGAGCATCCAAACCCTCAGACCTGCAAGGGG
CCAGCAGGGACGTGGCCCCACGCCACACACAACCTCTCCACATGCCTCAGCGCTGTTACTTGAATGCCT
TCCTTGAGGGAAGAGGCCCTTGAGTCACAGACCCACAGACGTGAGGGCCAGGGAGAGACCTAGGGGGTC
CCCTGGCCTGGATCCCCATGGTATGCTTGAATCTGCTCCCTGAACTTCTGCCAGTGCCTCCCCGTACC
CCAAAACAATGTACCATGGTTACCACCTACCCAGAAGACTGTTCCCTCCTCCAAAGACCCTTGTCTGC
AGTGGTGTCTCTGCAGGCTGCCCGTTAAGATGGTGGCGGCACACGCTCCCTCCCGCAGCACCACGCCAG
CTGGTGGCGCCCCACTCTGTCTTCTTCACTTACAGACAAAGGATTCTCAACCTTTGGTCAGTTA
ACTTGAAAACCTTTGATTTTCAGTGCAATGACTTTTTAAAGACACTATATTGGAGTCTCTTTCTCAGA
CTTCTCAGCGCAGGATGTAATAGCACTAACGATCGACTGGAACAAAGTGACCGCTGTGTAACACTAC
TGCCTTGCCACTCACTGTTGTATACATTTCTTATTTACGATTTTCATTTGTTATATATATATAAATA
TACTGTATATATATGCAACATTTTATATTTTTCATGGATATGTTTTTATCATTTCAAAAAATGTGTATT
TCACATTTCTTGGACTTTTTTTAGCTGTTATTCAGTGATGCATTTTGTATACTCACGTGGTATTTAGTA
ATAAAAACTATCTATGTATTACGTACATTA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_021180.4
Summary:	This gene encodes a member of the grainyhead family of transcription factors. The encoded protein may function as a transcription factor during development, and has been shown to stimulate migration of endothelial cells. Multiple transcript variants encoding distinct isoforms have been identified for this gene.[provided by RefSeq, Aug 2010]
Locus ID:	57822
MW:	33.1