

Product datasheet for **SC210406**

ARHGEF35 (NM_001003702) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ARHGEF35

Synonyms: ARHGEF5L

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_001003702

Insert Size: 876 bp

Insert Sequence: >SC210406 3'UTR clone of NM_001003702
The sequence shown below is from the reference sequence of NM_001003702. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GACAAGGAAAACTTCTATCAGTACTTATATGTCACAGTTTCATGTTTATCCAGTTCAATGTATTTT
TAAATTTTCTTGAGACTTCTTTGACTGATAGATTATTGTGAAGTGTGTTTTAAATTTCCAAATGTT
TAGGGATTTTCATATCTTTCTATGCTGATTTCGAATTGGATTCCCTACAATGATTTCTGGTTTTCATC
TGCTCTGGATGATTACTATCTCTTTTAAATTTGTTGTGGCGAGTTTATAGGCCTAGGACAGCTCTATCT
TGCCATGTGTTTCGTCAGCACTCAAAAAAATATGTGTTTCTGCTGTTACTGTGTGGAATATTCTGTA
AATGCCAAATAGATTCTTTGGTTAATGGCATTGTTGAGTTGTTTATATTCTTGTGATTTTCTTTTTT
TTCTCTTTTTTAAAAACATGAAATGCTTCACAAAATTTGTGTAATCTATATTCAAGAACCATTTTAAT
CTTGTTGGCATCATTCTGATTTTAGAAAGTGCTGCCAAAGCCAACACTCCTTGCTGATTTTCTGTCTAG
TTCTATCAATTATTGAAAAAGGCATGAGGAAATTTCTAACTATAATTGTGGATTCATCTATATCTCCTT
TCAGTTTTCTTTTTTAAAAAAATCAATGAAATTTATTAATTAAGCATAAAGTTACTTTACATTT
ATCTACAACCACAGTGAATACAGTTCTTGGCATGAAGACACCACAACCTTTAGAAATTTAAAGCCTCCTC
ACCTGCAAGATTACATATATAAACTCCCACTATTGTTTCTCTAAGAGTGGATTAGTTCACCAAGTTAA
AAGTTATATGATCTAGAATATAATAAATGGAATGATTTACTCATAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001003702.3</u>
Locus ID:	445328
MW:	34.2