

Product datasheet for **SC210417**

CCDC59 (NM_014167) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CCDC59
Synonyms:	BR22; HSPC128; TAP26
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_014167
Insert Size:	878 bp

Insert Sequence: >SC210417 3'UTR clone of NM_014167
The sequence shown below is from the reference sequence of NM_014167. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTTCTTCAAAAAATACAAGAAAAATGTAAACATTTTGTCTACAGGTTAAAAATATCTGCTGCCTATT
AGGTTCTTCTGTGACATGTGCCTCCAGCAGTGAACATAATTTGTGACATAAACTGGATTGCTAACT
ATGCTAAATATAAGATGTTACATATTTTATTATGGTAAAAATTTCTAAATATGTTCTACATGTTT
CTTATTTATTTGCCTCTGAAGGAAGGTTGGCCTGAAGAACTGAAAGAACCTCTTATTTTGCAAGACAGG
CCCAAGCATGTAATACTTTGTACCATATGAGATTTATATGAAATAAATTTTAAAAATAAGGAATCA
GAGCTATCAATGAAGCATTTCAATGAAATATTTCAATTTAGTAAACAGTTACGTTGTTTTAAATTTAT
TTTAATGATGAGGGAGGCAAGACTCCTCTGGACTTTTTATTTTAAAGTACATGTTACAGGGTGA
TGGGTGCTGAAATAAACTGCCCTCACCATAAACAAGCGACCTTCACAAAGCAGTGAACCTCAAAGTCT
GATATTTATAAACATTCCCTTAAATGGATGTACATGGCTTGTAGTCTAAGGCTTTTAAATAATGATG
GAAACCCACAGTGTACAGCAAGAATTAAGCTTTCAGTTCCTATGGACCCTATACTTTAAGAACATT
ATCTCTGTCTCCCTTGTGGTACCCTTGTGAGCTGAGATAAGCCACTGGGATAACTGCTTCTCATGCTT
TGCAAATTGCAGAAGATGAGAAAGGCTTGAAATATATATAGTCATTGTAAGCAAAAAGGGCTGACC
AATAATACTTGCCAAATATTGAATAATAAATTTTATAAATGATAGCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_014167.5</u>
Summary:	Component of the transcription complexes of the pulmonary surfactant-associated protein-B (SFTPB) and -C (SFTPC). Enhances homeobox protein Nkx-2.1-activated SFTPB and SFTPC promoter activities.[UniProtKB/Swiss-Prot Function]
Locus ID:	29080
MW:	34.7