

Product datasheet for **SC210420**

KIAA1370 (FAM214A) (NM_019600) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: KIAA1370

Synonyms: KIAA1370

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_019600

Insert Size: 863 bp

Insert Sequence: >SC210420 3'UTR clone of NM_019600
The sequence shown below is from the reference sequence of NM_019600. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ACAAACCCCTCAGTTTTACCAAGATGTGATAAGGAGTGATGATTTAAAGTATTTACTCAGTACCCAAG
TTTGCAAGTAAAAATTAGCATAGAATGGAGTGACCAAATTAACAATCAGGAGAGTGGATTCTCTCTG
TTATCCTGGACAGTTTTATGAAAGGATTCTGAAATGAAATCCATATATTCATGTAGACTGGAAAA
ACTCATGTCCTAATCCTTTTTGTACTGTTGAAACCACTTCATTGGACATGTTGCAATAGCAAAACCCCC
AGTTAGATTAGTGTACACATTTCTCAGTTATTTAATATTTAATGTTTTCTTAATACTCAAGTGAT
GTTTGTCTCTAGTGTCTAATGTAGCACAAATCCTATGTAATCATACTATGTATTTTGACATTAAT
GTTGAAATCAAATATATGCACAAGTCTTTAATTTGTGTAATGTGTTAAGTGCTGTTTCAATTAAGTTAT
TGAAATGAGAATAAATGTTGAGCTTCTTTAAAGTAACACACTATGCAAGCATGTGACTTTTTATA
TCTCTCATGTTTAGTTTTATAACACCATATCCAGGTTGCTATCTCACATAGTAGTCCTTTAACATATT
GTATTAGCAGTGCAATGTGGACTAAGCTGCTTCACTTTCCCTTTGCAAGTTCAGATCATCATGCCATT
CATAGCCAGGATTCCTTATCCCCAAACAGTTCTATTTTCTTAATCACTACTATAGAGTCTTTACAT
TAAATTACTGTCGTATGCTAGATAATTTCTCAAATTGTTAAAGAATATGTACTTTGGAACAAATTA
GTATTTATATTGTAATATATTCAAAAAAACTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_019600.4</u>
Locus ID:	56204
MW:	33.1