

Product datasheet for **SC210423**

C2orf65 (M1AP) (NM_138804) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C2orf65 (M1AP) (NM_138804) Human 3' UTR Clone
Symbol:	C2orf65
Synonyms:	C2orf65; D6Mm5e; SPATA37; SPGF48
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_138804
Insert Size:	854 bp
Insert Sequence:	>SC210423 3' UTR clone of NM_138804 The sequence shown below is from the reference sequence of NM_138804. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

TTCTTCCTGCCTTCAGAGTGGGAGAAGGATCCCTCAAGGCCCT**TAAG**TCACCAGCACCAGAGCCCAGCTGC
 CCAGCTTAACCATATCCATGCTCAGGTTACATAATGGCTATCTGTGGTCAGACTTGCTCTATCCGCC
 TGAGCCTCTGTGAGTGAGGGCTGACTGGGAAACAACAGCCTTCCTGTCCTGTTTCAGTGCTGTCCCACTC
 CTCAAGTCTGGAAGCGACACACCCGAGCCTGTCTTTCTCCAGCAAGGACTTTTCATTTCTTTAGAAATCA
 TTTGCTACTGTTACACAGGTGAAGATTAACACCCAGTAAGCTTCTACCATTGTTAGGAGCATTTCATAA
 CTCAGAATTTCTTCTGTAGCTCTGTGAAGCAGGTGGATGAGGTGAGTACACCTTTGGTAAACTGGACC
 TCAGGAACAAGGATGAGGTTTTGAAAGCTCATAAAGACAAGTAAGATTGAAATCCAAGCCTCATTTTCAG
 AGCCTGTGCCCTTCCCACTACACCACAGGCTTCAGCCTCCAAAGAGACAAGTGCTTGGTACCTACATGC
 AAAGTGTGTGTGCTGGGGGTGGGAGGGCTGCCAGAACAGGGGAGAGGATGGTGTAAGAAAAAGACCTAC
 TCCTTTCTGTTACCTCTCTCCACATGTACCAACCTTCCTGTTGCTCCCTCCATCCACAGAATAATAGC
 TACCATTATATAAATGTTTACTCTGGGCTGGGAGCAGTGGCTCACACCTGTAATCCCAACACTTTGAGAG
 GCTGAGGTGGGATGATCACTTGAGGCCAGGAGTTTCAGACACAGCCTGAGCAACACTGTGAGACCCCCCGG
 CCATCTCTACATAA

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_138804.3
Summary:	This gene encodes a protein that is likely to function in progression of meiosis. A similar protein in mouse plays a role in gametogenesis in both sexes. Alternate splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2013]
Locus ID:	130951