

Product datasheet for **SC210427**

SDAD1 (NM_018115) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	SDAD1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_018115
Insert Size:	865 bp
Insert Sequence:	>SC210427 3'UTR clone of NM_018115 The sequence shown below is from the reference sequence of NM_018115. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTTTTGAAAAAGAGAAAAAGAATGAAGTAACCTTCCTGGCAAGTTTTCCATTCTAGAGAAGTCTAAGT
TTGTGTCCTTGCTCTGAAAATTGGTAAATCAAGCATGTTTGTTTACATTAAGTCCAGACACACTGT
ATTGTGAAAAGTCTGAACATGTGGCAGCAATTTGTGTTTTTATTTTGGAGACGGCTAATGGTAGGAA
TGTTAATGTAATAGTGGTGGTAAATGTAATCATTTTATTATCATTCATGCAAAAAAGTATGTAT
TGAGTGCCTATTTCATTGCTACTGTAGATGCAAAACGAATGAGCTGTAACCCCTTCACTCAAGGCATTGA
CAGCTTAGCTGTGAGGGTGGACACACATATGTGTATTTACAGTGCAGTGTAATAGTTCTGTAGTAGAG
GTAACTCCATATTTCTGTGAGGTCACTGAGGCCTTATGGACTAACTCTGTGAGGATAGGAGTTATATA
TTCTTATAAGACAAAACAAAACAGGACAATGTTACAAGAGTAAGAGGTTCTTACTTGTACATAGGCTTT
CCTGCTGAAAACAGGCCCTGCTGTACAGATTTTGGGTACATAATTTAGCTCTTTTAGTCAATCCAAGA
GATTTAAGTGACCCCCCCCCCGTGTGTTTTTGTGTTTTGTTTTGATGCAATGCAATGCAATGCAATGCA
TTTTTGGTTAAGACCTCACTTTTAAACTGCCTTAAGTATAAATAGTACCTTTGGAATATATTTAGTTC
ATCATTTGAGCTGCCTTCATACTGGTTTCTCAGCCTTCCTTCAGCCTGTAATATTTTCAGCCCACTGT
TTACCTTGTCTCAATAAAGGTTTCTAATGCCAATA
ACGCGTAAGCGGCCGCGCATCTAGATTTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018115.4</u>
Summary:	Required for 60S pre-ribosomal subunits export to the cytoplasm.[UniProtKB/Swiss-Prot Function]
Locus ID:	55153
MW:	33.4