

Product datasheet for **SC210438**

RAP1 (TERF2IP) (NM_018975) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RAP1 (TERF2IP) (NM_018975) Human 3' UTR Clone
Symbol:	RAP1
Synonyms:	DRIP5; RAP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018975
Insert Size:	851 bp
Insert Sequence:	>SC210438 3'UTR clone of NM_018975 The sequence shown below is from the reference sequence of NM_018975. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCTCGGAGGATTGAATTCGAAAGAAATATTGGCAAGATAATGAGAAAAGAAAAAGTCATGGTAGGT
GAGGTGGTTAAAAAAATTGTGACCAATGAACCTTAGAGAGTTCTGCATTGGAAGTGGCACTTATTTT
CTGACCATCGCTGCTGTGCTCTGTGAGTCCTAGATTTTTGTAGCCAAGCAGAGTTGTAGAGGGGATA
AAAAGAAAAGAAATTGGATGTATTTACAGCTGTCCTTGAACAAGTATCAATGTGTTTATGAAAGGAAGA
TCTAAATCAGACAGGAGTTGGTCTACATAGTAGTAATCCATTGTTGGAATGGAACCTTGTATAGTAG
TGACAAAGTGAAAGGAAATTTAGGAGGCATAGGCCATTTAGGCAGCATAAGTAATCTCCTGTCCTTTG
GCAGAAGCTCCTTTAGATTGGGATAGATTCCAAATAAAGAATCTAGAAATAGGAGAAGATTTAATTATG
AGGCCTTGAACACGGATTATCCCCAAACCCTTGTCATTTCCCCCAGTGAGCTCTGATTTCTAGACTGCT
TTGAAATGCTGTATTTCATTTTGCTAAGTATTTGGGTACCTGCTCTTTGGCTGTTCTTTTTTTG
GAGCCCTTCTCAGTCAAGTCTGCCGGATGCTTTTCTTACCTACCCCTCAGTTTTCTTAAACGCGCA
CACAACCTAGAGAGTGTTAAGAATAATGTTACTTGGTTAATGTGTTATTTATTGAGTATTGTTTGTGC
TAAGCATTGTGTAGATTTAAAAAATTAGTGGATTGACTCCACTTGTGTGTTGTTTCATTGTTGAA
AATAAATATAACTTTGTATTCGA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018975.4</u>
Summary:	The gene encodes a protein that is part of a complex involved in telomere length regulation. Pseudogenes are present on chromosomes 5 and 22. [provided by RefSeq, Apr 2010]
Locus ID:	54386
MW:	32.2