

Product datasheet for **SC210439**

ZNF676 (NM_001001411) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF676
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001411
Insert Size:	850 bp
Insert Sequence:	>SC210439 3' UTR clone of NM_001001411 The sequence shown below is from the reference sequence of NM_001001411. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

TTCATACTGGAGAGAATCCCTAAAAATGTGAAGAATGTGGTAAATCCTTCAGCTGATTCTCAGTCTAGTA
AACATAAGAAAAATTCATACTGGAGAGAAAAATCTACAAAGGTGAAGAATGTGGCAAAGCTTTTAAACAGTC
CTCAATCTTTATTACACATAAGAAAAATTCATGCTGGAGAGAAACACTACAAACGTGAAGAATGTGACAAA
GGCTTCAACCGGTCCTCAAACCTTGTGGAACATAAGAGAATTCATACTGGAGAGAAACCTACAAATGCG
AAGAATGTGATAAAGACTTCAACTGTTCTCACACCTTACTACTCATAAAGAATTCTACTGGAGGAAA
AGCCCTACAAATGTGAAAAACGTGGCAAAGCTTTTCATTAGTTCTCAACCTTACTGAACATAAGGGAAT
TTATACAGGAGGGAACCTACAAATGTGAAGAATGTGGCAAGCCTTTTAGCTGTTCTCAATCCTTACT
AAACATTTCATGTAATTCATACTGGAGGGAATTCCTATAATTGTGTGAATGTGGCAAAGCCTTTAACCAG
TCCTTAAGGCTTACTACATATAAGACAACCTACTGGAGAGAAACCATGCATGTGTGAAGAATGTGGGA
AAGCCTCTAACAGATCCTCAATTCTTAAGAAACATAGGCTAATTCATACATGAGAGACTCTACAAACCTG
AAAGATGTGACAATGCTTTTGGCAACACCTCAGACTTTTCTGAATATAAAGAAATCATATGAATGAGAA
ATCCTAGAAATGTGAAGAATGTGACAAAGCCTTTAAATGGCTGTCACACTTGATTGTAGGTAAGATAATA
TACTGGAGAA

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001001411.2</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	163223