

Product datasheet for **SC210444**

RFESD (NM_001131065) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	RFESD
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001131065
Insert Size:	853 bp
Insert Sequence:	>SC210444 3' UTR clone of NM_001131065 The sequence shown below is from the reference sequence of NM_001131065. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

TGCCACTGGAGACTTCAAAGTAATTAAGAGTTCTT**CTCGA**TAAAAATATATAGAAATGAAAAATGTTGT
GTATGCTTGAAAACATTTTGAATAACTCTGCTTCAGTTTTAAAGGTGATGAAATTTTCAACATAGGC
ACAACGTGTTTAAATTCACATACATTTGAGAGGTTAAGAGCACACATACTTAGTATGATGTAAGGATTA
TCATCAGGATCAATAGAATACATTTATCGAATCTTCTGGATTAATTAGAACTGAAGGTTATAGGTTT
GGGATGTTCTGTTAAGAAGGATAAAAAATAATTTGGAGTAAACAAATAATAGGTAATAAAGGATTA
GTAAGAAAAAATTTTAAATTTAAAAAATTTAATTTTAAAAATAGGAAAAATCCCTGTTTGATGTAGCA
GACTTAATCTGGAAAAGTTTTAGTATATACAAAGTAAACAGACTATATAATGAACTTCCATCACTCAGC
TTTAACAGTTATCAGCATTCTCCATTCTGTTTTATTATACCTCTGTGGATTGCTACCTCTCCTATT
TGATTGTTATTTTACATTTTATAGGAAAGTTTAAATACATTGAAATGTTACCTCTGAACTGTGCAACCT
TCCTAAATAACAGACTTTTATACTTAAAAATGCCTCCAAAAATGAGATAGGTTTACCTGCTGAACTTG
AAAATAAGAACTTAACCACTGTTAAGTGATCTGTGATCAAGGCAGATCTTGAAGTAAGGAAGCCAGACA
CAGCTGCTGCCAGCCATGAGAGATGTAAGTCTCCTATTCAGTGAGCCAATGTTCTTTGAATGCCTATA
TGAAATTTTCAAA

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001131065.1</u>
Locus ID:	317671