

Product datasheet for **SC210445**

RFESD (NM_173362) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: RFESD

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PS100062)

ACCN: NM_173362

Insert Size: 853 bp

Insert Sequence: >SC210445 3' UTR clone of NM_173362
The sequence shown below is from the reference sequence of NM_173362. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site **Blue**=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

TGCCACTGGAGACTTCAAAGTAATTAAGAGTTCTTCC**TGA**TAAAAATATATAGAAATGAAAAATGTTGT
GTATGCTTGAAAACATTTTGAATAACTCTGCTTCAGTTTTAAAGGTGATGAAATTTTCAACATAGGC
ACAACGTGTTTAAATTCACATACATTTGAGAGGTTAAGAGCACACATACTTAGTATGATGTAAGGATTA
TCATCAGGATCAATAGAATACATTTATCGAATCTTCTGGATTAATTAGAACTGAAGGTTATAGGTTT
GGGATGTTCTGTTAAGAAGGATAAAAAATAATTTGGAGTAAACAAATAATAGGTAATAAAGGATTA
GTAAGAAAAAATTTTAAATTTAAAAAATTTAATTTTAAAAATAGGAAAAATCCCTGTTTGATGTAGCA
GACTTAATTTCTGGAAGTTTTAGTATATACAAAGTAAACAGACTATATAATGAACCTCCATCACTCAGC
TTTAACAGTTATCAGCATTCTCCATTCTGTTTTATTTATACCTCTGTGGATTGCTACCTCTCCTATT
TGATTGTTATTTTACATTTTATAGGAAAGTTAATAACATTGAAATGTTACCTCTGAACTGTGCAACCT
TCCTAAATAACAGACTTTTATACTTAAAAATGCCTCCAAAAATGAGATAGGTTTACCTGCTGAACTTG
AAAATAAGAACTTAACCACTGTTAAGTGATCTGTGATCAAGGCAGATCTTGAAGTAAGGAAGCCAGACA
CAGCTGCTGCCAGCCATGAGAGATGTAAAGTCTCCTATTCAGTGAGCCAATGTTCTTTGAATGCCTATA
TGAAATTTTCAAA

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_173362.3</u>
Locus ID:	317671