

Product datasheet for **SC210449**

RTL8A (NM_001078172) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RTL8A (NM_001078172) Human 3' UTR Clone
Symbol:	RTL8A
Synonyms:	CXX1b; FAM127B; MAR8A; SIRH6
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001078172
Insert Size:	865 bp
Insert Sequence:	>SC210449 3'UTR clone of NM_001078172 The sequence shown below is from the reference sequence of NM_001078172. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TTTGGATGGGAGGAGGACGAGGACTTCTAGGCCGGGAGACCCTTGGGCCTGGGGCGGGTGCTCTGGGA
AGAGTTCGCTGTGCCAGTGGCCACCGCTAGGGTCTCCACAGGCGCCCTCCCTCCGCGCTCCCTCCCCC
TCGAGCCGCCGCGATGTCCCTGCGCTCCTGTCCCTCCCGCTAGTGCTTGCCTTTGTTCCAGGAATA
GCGCTCCAGGCTCCTGCTGCAGCTCCTGGGCCGCACTCTAGAGCGCGCTGCCTCCGCCCTCTTTTGCCA
GCCCAGCCCTCCATGCACATTTGGACGCTGTCTGCAGTTCAGGTGCAAGCTGGGCTCCTGTTACATA
CTGGACAGACCACCCACAGCCGCTGCTGCCAACCCACCTCCACTCCTCCAGACTGCCAGCCAACATACAT
CATTCTGCCCACAGACCTACGCTACCACTGCCATCGCCATCCACCGCATCTCACCAACAGACTGTGGCT
CCTAGTGATGTGGACTCACCTCGGAGGTATCCGAGCTGGACACAGCCCTGGACAGTGATCCAGACAGC
TGGCCGTCCCCAAGGAGCCTGTACCTTCAGCGAGACCCATTTCTCCCATCCCCAGAGACCTCTTG
TGTTCTGCCACATAGCTGCCAGGGCTTAAGTGTGCCTGGCAACCAATCGAATCTCTCATTCTCTCCT
GTGGACCAGTTAGTTTTGCCTAGAATCCTGTTTCTTGAATTTGCAGCTGTCTCTGATGGGTGCC
TTTTGTTCAACACAGTAAGCCCTGCTCCCTTCCCTGCTCTAATACTACCTGTACAAAGGTTTTTCC
TTATTTTAAATAAATGTAGACACTATTAATAGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001078172.2</u>
Locus ID:	26071
MW:	31.1