

Product datasheet for **SC210459**

ZNF155 (NM_198089) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF155

Synonyms: pHZ-96

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_198089

Insert Size: 882 bp

Insert Sequence: >SC210459 3'UTR clone of NM_198089
The sequence shown below is from the reference sequence of NM_198089. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CTTTTATCATTATTTCTAAATGACACATAAATTGTTTCATATTTATGGGGTACAACGTGCTATTTTAATGT
GTGCATACAATTTATAGTGATCCAATCAGTGTAATTGGTGTATCTGTTACCTCAAACATTTACCATTTT
TTTGTGTTGGAAAATTCAAAATCCATTTGAGAGGAGTTGGTAGACAGCAAGGGAAGGGTCCCTGGAGAC
CCCCAGCCAAACGGGTCAAGTGTCTCATCCCCACATAACACAAAAAGCAGCCTGGGGAAGACAATCAAGC
TACAGGCACTGATAAGGGAAGTACACAGGGGATTGTGCTGGAGACATGCCCATGGCTGCACAGATAG
AAGAACCTCCAGCCATTCAAGTAAAACTCGCACAAACCTCCGGCTCACTCAGTTAAGGGAACAAGAC
CAGACACAGAAATGCCTTTGTCCTTTGTATAATCAGCAGGCTCCGAGGAAAAAGTTTCTTCTTTTGTGG
GCATGAACACAGTGGGCTTCAGTAGGTTCAAGTTGGCACTTTCTTTTGTGAAAAGGAAAAGCCAGCC
TCTGTAATCATCACTTCAGCTCCTGATTGGTCCCAGGCCAAGCTGAATAGCCCTTAGGAATAATCACT
TCAACTCCTGATTGGTCCCGGCCAAGCTGAATAGCCCTTAGGAATCATCACTTCGGTTCCTGATTGGT
CCCGGCCAAGGTCCAGACCAAGCTGAATCAGCCTTTCTCCAAGACAGCCCGCAGGCTAAGTGCAATTC
CTTCCCTTCCAGTCCATAAAAAACCCGACCCAGCCTCATACTGGGCAACTATTGGATCCCCAT
CTCCACTGTGGAGAGCTTTTTCTTCACTTATTAACCTTTTACTCCAATTTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_198089.3</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	7711
MW:	33.5