

Product datasheet for SC210460

ZNF155 (NM_003445) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF155 (NM_003445) Human 3' UTR Clone
Symbol:	ZNF155
Synonyms:	pHZ-96
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003445
Insert Size:	882 bp
Insert Sequence:	>SC210460 3'UTR clone of NM_003445 The sequence shown below is from the reference sequence of NM_003445. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTTTTATCATTATTTCTAAATGACACAATTGTTTCATATTTATGGGGTACAACGTGCTATTTTAATGT
GTGCATACAATTTATAGTGATCCAATCAGTGAATTGGTGTATCTGTTACCTCAAACATTTACCATTTC
TTTGTGTTGGAAAATTCAAAATCCATTTGAGAGGAGTTGGTAGACAGCAAGGGAAGGGTCCCTGGAGAC
CCCCAGCCAAACGGGTCAAGTGTCTCATCCACATAACACAAAAAGCAGCCTGGGGAAGACAATCAAGC
TACAGGCACTGATAAGGGAAGTACACAGGGGATTGTGCCTGGAGACATGCCCATGGCTGCACAGATAG
AAGAACCTCCAGCCATTCAGGTAACAACTCGCACAAACCTCCGGCTCACTCAGTTAAGGGAACAAGAC
CAGACACAGAAATGCCTTTGTCTTTGTATAATCAGCAGGCTCCGAGGAAAAAGTTTCTTCTTTGTGG
GCATGAACACAGTGGGCTTCAGTAGGTTCAAGTTGGCACTTTCTTTTGTGAAAAGGAAAAGCCAGCC
TCTGTAATCATCACTTCAGTCTCTGATTGGTCCCAGGCCAAGCTGAATAGCCCTTAGGAATAATCACT
TCAACTCCTGATTGGTCCCGGCCAAGCTGAATAGCCCTTATGAATCATCACTTCGGTTCCTGATTGGT
CCCGGGCCAAGTCCCAGACCAAGCTGAATCACGCTTTCTCCAAGACAGCCCGCAGGCTAAGTGCATTTC
CTTCCCCTTCCAGTCCATAAAACCCCGACCCAGCCTCATACTGGGCAACTTATTTGGATCCCCAT
CTCCACTGTGGAGAGCTTTTTCTTTCACTTATTAACTTTTACTCCAATTTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_003445.4</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	7711
MW:	33.5