

Product datasheet for **SC210472**

COPS7A (NM_001164093) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: COPS7A (NM_001164093) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: COPS7A
Synonyms: CSN7; CSN7A; SGN7a
ACCN: NM_001164093
Insert Size: 839 bp
Insert Sequence: >SC210472 3' UTR clone of NM_001164093
The sequence shown below is from the reference sequence of NM_001164093. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

CGCCAAGATTTGGTCCAAGTCGAAT**TGA**AAGGACTGTCGTTTCCTCCCTGGGGATGTGGGGTCCCAGCTG
CCTGCCTGCCTCTTAGGAGTCCCTCAGAGAGCCTTCTGTGCCCTGGCCAGCTGATAATCCTAGGTTTCATG
ACCCCTTACCTCCCCTAACCCCAACATAGATCACACCTTCTCTAGGGAGGAGGCAAATGTAGGTCATGT
TTTTGTTGGTACTTTCTGTTTTTGTGACTTCATGTGTTCCATTGCTCCCGCTGCCATGCTCTCTCCCT
TGTTTCCTTAAGAGCTCAGCATCTGTCCCTGTTTATTACATGTCATTGAGTAGGTGGGTAGCCCTGATGG
GGGTCGCTCTGTCTGGAGCATAACCCACAGGCGTTTTTCTGCCACCCATCCCTGCATGCCTGATCCCC
AGTTCCTATACCCTACCCTGACCTATTGAGCAGCCTCTGAAGAGCCATAGGGCCCCACCTTTACTCAC
ACCCTGAGAATTCTGGGAGCCAGTCTGCCATGCCAGGAGTCACTGGACATGTTTCATCCTAGAATCCTGTC
ACACTACAGTCATTTCTTTCTCTCTGTCCTTGGGTCCTGGGAATGCTGCTGCTTCAACCCAGAG
CCTAAGAATGGCAGCCGTTTCTAACATGTTGAGAGATGATTCTTTCTTGGCCCTGGCCATCTCGGGAAG
CTTGATGGCAATCCTGGAAGGGTTAATCTCCTTTTGTGAGTTTGGTGGGGAAGGGAAGGGTATATAGAT
TGTATTAATAAAAAAAAAAGGTATATATGCATATATCTATATAATAATATGACGCAGAAATAAATCTATGAG

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001164093.1
Summary:	This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described. [provided by RefSeq, Mar 2014]
Locus ID:	50813