

Product datasheet for **SC210480**

HID1 (NM_030630) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HID1 (NM_030630) Human 3' UTR Clone
Symbol:	HID1
Synonyms:	17orf28; C17orf28; DMC1; HID-1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_030630
Insert Size:	858 bp
Insert Sequence:	>SC210480 3'UTR clone of NM_030630 The sequence shown below is from the reference sequence of NM_030630. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTGAAGCTGTTTGAGATACAGCGGGTGTGAGGATGAAGCCGACGAGGGGCTCAGTCTAGGGGAAGGCAG
GGCCTTGGTCCCTGAGGCTTCCCCATCCACCATTCTGAGCTTTAAATTACCACGATCAGGGCCTGGAA
CAGGCAGAGTGGCCTGAGTGTCATGCCCTAGAGACCCCTGTGGCCAGGACAATGTGAATGGCTCAGA
TCCCCCTCAACCCCTAGGCTGGACTCACAGGAGCCCCATCTCTGGGGCTATGCCCCACCCAGAGACCAC
TGCCCCAACACTCGGACTCCCTCTTTAAGACCTGGCTCAGTGTCTGGCCCTCAGTGCCCACTCC
TGTGCTACCCAGCCCCAGAGGCAGAAGCCAATGGGTCACTGTGCCCTAAGGGGTTTGACCAGGGAACCA
CGGGCTGTCCCTTGAGGTGCCTGGACAGGGTAAGGGGGTGCTTCCAGCCTCCTAACCCAAAGCCAGCTG
TTCCAGGCTCCAGGGGAAAAAGGTGTGGCCAGGCTGCTCCTCGAGGAGGCTGGGAGCTGGCCGACTGCA
AAAGCCAGACTGGGGCACCTCCCGTATCCTTGGGGCATGGTGTGGGGTGGTGAGGGTCTCCTGCTATAT
TCTCCTGGATCCATGGAATAGCCTGGCTCCCTCTTACCCAGTAATGAGGGGCAGGGAAGGGAAGTGGG
AGGCAGCCGTTTAGTCTCCCTGCCCTGCCACTGCCTGGATGGGGCGATGCCACCCCTCATCCTTAC
CCAGCTCTGGCTCTGGGTCCCACCCAGCCCCCGTGTGAGAACAATCTTTGCTCTGTACAATCGG
CCTCTTTACAATAAAACCTCCTGCTCCACA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCACTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_030630.3
Summary:	May play an important role in the development of cancers in a broad range of tissues. [UniProtKB/Swiss-Prot Function]
Locus ID:	283987
MW:	31.1