

Product datasheet for **SC210541**

WDR21A (DCAF4) (NM_001163509) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: WDR21A

Synonyms: WDR21; WDR21A

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_001163509

Insert Size: 928 bp

Insert Sequence: >SC210541 3'UTR clone of NM_001163509
The sequence shown below is from the reference sequence of NM_001163509. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CAGGACCTTTACTGTTACTCCTACAGCTAATTCTGCAGGGCACAGCCAGAGCCATGTGGATTTGACTT
ACGGGAGTAAAGCGTAACCTTTTACTGCATCTAATGAGGGTGTCTTAAGTGACACTCAGTGTACACAGA
TCCCATCCTCTGGCTGCTAGGAGAGAAGTGCTGAATGTTCCGTGTGGAGATGCTCAGGAAAGTTATTTG
AGTTAAATTGCTGGCTGAGAGAGCTTGGAAGTCCTTTTCATAAAAGGTACCTCTTTCCTTTTCTTATTG
AATTCTTAGAACTTAGTTAACCTCCCTGCCTTTTCTTAACAAAAGGACTTTTCTAAGGACTGAAGAT
TGGCAAAACGAAAGCTTCTCCTCCAAGAGCCATTGAAGAAGCCAGTGATGAGACGGTGAGATGG
TTTGAGTCCTCGGTGCCTGGGTAGCAGGAAGAAAGACCTGCATCCTGCATCTGTAATTGGGAAGCCAG
CGGAGAGGACGGGGAGGTTACTTCTAAGTTTCTGCAGAAATATTGAAGGCTGGAGTTTGGAAATCCTT
AAACTTGGCCTTCTCAAACCTCAGCAGCAGATCTCCGGGATTCTGCTGTTATTATCCAAGGCGTTGGAA
GGAAAGATGGATCTTCTACATGCTAGAAGTTTAAACGGTCCTTAACATGCCTTTGTTCAAGCACCTT
CCAGAAATGAAGGTTACGACGCTCTGGTTTCTATTACGGTGACTTGAATGTCAGATTCAAGGGCCCGGC
GTCAAAGGAAATTGGTTTTGACTTTTTGTAATCTAGGAGCGACAGTTTCGTGAGATGTTTATTCAGTGTT
AAAGAGCCTGTTTTCTACCAACAATAAAACCAAGAGAAGAATCATGGTGTGACATGCACTTCTGTGA
ATGACTTAGGAGATCATGAGCTTGGGCTGCT
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001163509.2</u>
Summary:	This gene encodes a WD repeat-containing protein that interacts with the Cul4-Ddb1 E3 ligase macromolecular complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2009]
Locus ID:	26094
MW:	35.8