

Product datasheet for **SC210550**

FAM118B (NM_024556) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	FAM118B
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_024556
Insert Size:	865 bp

Insert Sequence: >SC210550 3'UTR clone of NM_024556
The sequence shown below is from the reference sequence of NM_024556. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CACAGTGAAATAAGAGGCTGTAGTACATGACGAGCTAGAGAAATCACCACCGTTTAGACCAAGCTGTA
AGGCCCTACTACAGACAGTGTTTAACAAGTAACTTACAAGAACCACACAATCCCAGAAAAGTAACA
ATAGCCAGAGGTTGAAGGGCGGGGTAGAAGAGGGGGGAATGTTGCAGCGTAATCCTTCATACCACCTGG
TTCTTGATATTCTGCCGCTGTTCAAGTTCAAGAATAAAAGCGACAGCAGGACCCAAATGCAGCTCCCA
ACCCACTCCCAGGCTAGACATGCTTGTGTCCACACAGCACACCAATGTGATACTTCCACTGACCGGCT
GCAGCTCTGCATGAAGGACTCGGGGTCTGGATGCCATGGAATCACTGTGGCTCTTGTTGCAGTTTTGTA
CTCTATACTTGTTTTTCAATTAAGCTTAATGGCTTTTTTAAACATGACTTGAAGCTCTAGTTTTCTA
GATCTTTTACAGGTACAGTATTTTACATAACTAAGCTGTATTAAGCTTGTTCATTTACTTGCCAGG
ACCCTGGCTCTACTTTTAGAGTCATTGTAAGAACTCTAAGTGCATCAAGGTACTAATAAGCTTAATT
TTAATAACCCAAAGTTTAAAGTTCCGATCTTTCTCCTTGGGGTGGAGTGATCTCATTCTCAGGACAAC
CGTTTACTTACCTGATTCTCGGAGCATTATCAACTTCTGCTCTGTTGTCTGACCATACATATGCTCT
AGAAGTACAGTTAAGTGTGTTGTGGAATTTAGTTTTGAATCCTGTATAAAGAACTCCAAGACCTAAA
AGGCAATTTATTTATGAAATTTATTTTCTATATAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024556.4</u>
Summary:	May play a role in Cajal bodies formation.[UniProtKB/Swiss-Prot Function]
Locus ID:	79607
MW:	33.1