

## Product datasheet for **SC210553**

### FAM185A (NM\_001145269) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** FAM185A (NM\_001145269) Human 3' UTR Clone  
**Symbol:** FAM185A  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pMirTarget (PS100062)  
**ACCN:** NM\_001145269  
**Insert Size:** 853 bp  
**Insert Sequence:** >SC210553 3'UTR clone of NM\_001145269  
The sequence shown below is from the reference sequence of NM\_001145269. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGGTTTCAGTCCCTGAAACTGCAGGACTAAAAACTGATTTGAGTTGGATTTATGATTTTTAACAAATGATT
CGCAGATCTGTGACTACAAAAATGTAATCCCACCATTATATAAAGGTTGAAAACAACAAATTGAGAA
TGAATACTGGTGAAGTGTGGGAGGCTTTTCAAATTTGTGCTTTGCTATTGTTATTTCTAC
TGCTCTGTAACAAATTACCACAAATTTAGCAGCATAAAAATAACTCATTACAGCTATGTCTAGGTTT
TCCGCTCAGGGAATCACAAAGGCTAGAAATCAAGATTTTGCCAGGCTGAGTTCTCATCTGGAGTCTCTG
GGAAGAAATCTACTTCCAATTCATTAGATTATTGGGAAAATCCAGTTCCTTGTGACTGGATGACTGA
GGTCTCCATTTTCTTGTGGCTGTGAGCCAGCAATCACTCTAAGCTTCTAGAGGCCACCCACATGCCTT
CTCATGTTGCCCTCCATCTTCAAGCCAGCAGTGATGCGCTGAGTTCTTCTCATGCTTTGAATATCTGA
CTTCTTCTGCAGCCAGCTGGAAGAATATATACTTTAAAGGGCTCACTTGATTAGGTCAGGCCACCCA
GGATAACCTTTGTTTGGCATATAATGTAACATTCATAGGAGTGATGGCCATCACATTCAGTTTTAC
TCCCACTCAAGGAAAGGAGATTATACAAAGGCATCCATCATGTGGGACAATTAGCTTAGAATCTTTA
TAATGTAGCTCTCACATGACAGGTTTTCTACTTAGTACCTTTTATGGCATCAATGTAATTTAAGCTGT
AAATAAAAAACAATATGTGAAGCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_001145269.2</a></u>
<b>Locus ID:</b>	222234
<b>MW:</b>	32.1