

Product datasheet for **SC210586**

WBSCR16 (RCC1L) (NM_030798) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: WBSCR16

Synonyms: WBSCR16

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_030798

Insert Size: 879 bp

Insert Sequence: >SC210586 3'UTR clone of NM_030798
The sequence shown below is from the reference sequence of NM_030798. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATGGTGACCTGGCCAAGTCATTATCTAAACCTCCCTCACCTGCTTGGGCGGCCCCGTCCCGGGAACC
ACTGGCACTCCTTGGCAGAGGCCAGCGCGTGGCCAGCCCCCGGGGTTCTTGGATGGTGGTGGCGGAGG
ACCTGCGTGCAAGTGACGCTCTGTCTGAATCCCTTAGCGGGTACCTACCAGGAGGATCAGGGCAAG
GTCCCTCTCCAGCTGCAGGTGAGGCCTGCGGAACCTCAGCTTGGATGGCAGCCTTTGGTGGGCGCGTGTG
GCGCGCACGTCTCTGTTCTCTCAAGTAACATGCGACGGTGTCTGGTGTACGTCTCGCCTGAGAAGCC
CGTCTTAGGAAAGCTTAGCTTGAACACAGTGCTCGGGAGGTTTCTGCTCTGTCTGTCTGTCATGGCAGTCTCT
TGGTTTGTGTCTGGCCAAGGCCATGCGTGTGCCTCGGACCGAGCCCAGCTTAGCGAGGGAGTCAGGC
TGGCTTCGGCCCTCGGTTTTTCATTACAGGCCACCCTGCTCATGGCCCTTCTGGCCGCGCTGCCACACCGC
AAGCTCGCTGGGGGACACTAGAAGCACCGTGGCCTGGGATTCCATCTGGAGCTGTCCGCAGGCACCCAG
CCCCAGCCTCCACACGCTCACTGCCTGGCTTGGAAAAGTAAAGAAGCCCTCAGGAAGAGAATCGAG
GCCAAGTTCCTCTGCGCCGAGGGCCCCGAGCATATCCGCCAAGGCTCAGCTGCAGTGCCAGGCGGAGGA
GGAAGATCCAGAAATTGTGAACAATGTTTGAATTTAGTAGCGTGACTTGCCCTTTCCCTTTAAAAACATCT
TTTACAAATCTGTCTTGAATAAAGTCTATTTTCTGCCTTTTGGTTTTTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_030798.5</u>
Summary:	This gene encodes a protein containing regulator of chromosome condensation 1-like repeats. The encoded protein may function as a guanine nucleotide exchange factor. This gene is located in a region of chromosome 7 that is deleted in Williams-Beuren syndrome, a multisystem developmental disorder. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]
Locus ID:	81554
MW:	32.1