

Product datasheet for **SC210593**

splicing factor 1 (SF1) (NM_201998) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	splicing factor 1 (SF1) (NM_201998) Human 3' UTR Clone
Symbol:	splicing factor 1
Synonyms:	BBP; D11S636; MBBP; ZCHC25; ZFM1; ZNF162
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_201998
Insert Size:	878 bp
Insert Sequence:	>SC210593 3'UTR clone of NM_201998 The sequence shown below is from the reference sequence of NM_201998. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGAGTGAGGACTTCCGCGCCATTGGTACCCCTTCCAGGCAGACAGCCTCAGCAACGCCCTGGTGGGA
CAGGATGGTTCGGCAAAGCAGCCTGAGTTATTTTGTGGACGGAATCGGAACACGCTGGCTCCATATCG
TGAAATTTTATTAATTTTTTCTTTTCTTTGTTACTTCTTTATCTTTTCTTTCTTCAGACTCCGT
CCAAGGAGATGCTCTCCCGGTCTTCTGCTGCAATTTAGATTCTTTCCCTTCTCTCCAGTTCTCCTTC
CCTTACCAAGGAGAGGGGAGCAATGGTTTTGGGCAAGGGCTTTGGCCATTTCATGTCAAGCTGGTTGTG
GGTTTTTCAAGGTGCCATAGCCACCCCAAGATGTTTGTAAAGCGTGGGGTTTTTAACTCTCTGCC
ACCCTTGTCAAGGGAGTCTTGTAAGATTGCCGAGGGTAGGTTTCATCTCCAGGTTTCGGGATTCCCATCC
GTCCTGGCGATCCTGCCAGCAGTGGGTGGGCAGCCTGAGCTCCCTCGGGCTCGCCTGCCAGCCTGGAGT
TCTTCTGTGCTCCTTGATCACCTGAGCTGCCTCAGATTCCATTTGGTCTCTCCTTCTGGAAGGCTT
CCTTTTATGTTTTGTTTTAATCCCAATGTCTGAATGTTTGCAGTGTGTAGGGGTTTGAGCCCTTGT
TCATTCTCCTTCTTTTCTCCCGCTTCCCTCTCCATGAAGTGATTCTGTTGACAATAATGTATACTG
CGGTTTCTCTTCACTGGTTTATCTGCAGAAATTTCTTGGGCTTTTTTCGGTGTAGATTCAACACTGC
GCTAAAGCGGGGATGTTCCATTGAATAAAGAGCAGTGTGGTTTTCTGGG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCGCCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_201998.3</u>
Summary:	This gene encodes a nuclear pre-mRNA splicing factor. The encoded protein specifically recognizes the intron branch point sequence at the 3' splice site, together with the large subunit of U2 auxiliary factor (U2AF), and is required for the early stages of spliceosome assembly. It also plays a role in nuclear pre-mRNA retention and transcriptional repression. The encoded protein contains an N-terminal U2AF ligand motif, a central hnRNP K homology motif and quaking 2 region which bind a key branch-site adenosine within the branch point sequence, a zinc knuckles domain, and a C-terminal proline-rich domain. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2016]
Locus ID:	7536
MW:	32.5