

## Product datasheet for **SC210604**

### **PDE8B (NM\_001029851) Human 3' UTR Clone**

#### **Product data:**

Product Type:	3' UTR Clones
Product Name:	PDE8B (NM_001029851) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	PDE8B
Synonyms:	ADSD; PPNAD3
ACCN:	NM_001029851
Insert Size:	1732 bp



[View online »](#)

**Insert Sequence:** >SC210604 3'UTR clone of NM\_001029851  
 The sequence shown below is from the reference sequence of NM\_001029851. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAAGTTTGAGGCTTCCATCTGACAGCTAAAGCCAAGCCACAGAGGGGCCTCTTGACCGACAAAGGAC
ACTGTGAATCACAGTAGCGTAAACGAGAGGCCTTCTTTCTAATGACAATGACAGGTATTGGTGAAGGA
GCTAATGTTAATATTTGACCTTGAATCATTCAAGTCCCAAAATTTATTCTTAGAAAAGTTATGTTCCA
TGAAGAAAAATATATGTTCTTTTGAATACTTAATGACAGAACAAATACTTGGCAAACCTTTTGCTCTG
CTGTATCCTGTGTACCCTTGTCAATCCATGGAGCTGGTTCAGTAACTAGCAGGCCACAGGAAGCAA
AGCCTTGGTGCCTGTGAGCTCATCTCCAGGATGGTGACTAAGTAGCTTAGCTAGTGATCAGCTCATCC
TTTACCATAAAAAGTCATCATTGCTGTTTAGCTTACTGTTTTCTCAAGAATCATGATCTGAAGGATTC
ATAAGGAGCTTATCTGAACAGATTTATCTAAGAAAAAAGGACATAAAATAAGTGAACAACATA
GGACCAAATTACAGATAAACTAGTTAGCTTACAGCCTCTATGGCTACATGGTCTTCTGGCCGATGGT
ATGACACCTAAGTTAGAACACAGCCTTGGCTGGTGGGTGCCCTCTCTAGACTGGTATCAGCAGCCTGTG
TAACCCCTTTCTGTAAAAGGGTTCATCTTAACAAAGTCATCCATGATGAGGGAAAAAGTGGCATTTC
ATTTTTGGGGAATCCATGAGCTTCTTTATTTCTGGCTCACAGAGGCAGCCACGAGGCACTACACCAAG
TATTATATAAAAGCCATTAATTTGAATGCCCTTGGACAAGCTTTTCTTAAAAAAGTT
TATATACATGTTAAAAATTTTATTAATAATCCAAATTTCTGGGTGTAAGCCATGCAGTATGTTGTGC
CATGCACCATCTATACGTAATATTTGGGAGGGGATAGGGGAAGTTTCAAGGTTTCAAGATTTTTAAACC
AGTCTATATTTGAAATCATGCCATTAGCTTTAACAATGTTAAATTTGTGGTGAATATATACAGAAA
ATAGGCATATGGGAGGCAGTAAATACTATTTTAAAGCTATTAATTGTATGCTAAAAGCTTCTAAAGCAC
AAGTGCATATTTTATACTGCTCTTTTCAAACTCCCTGTGATCTACATAAAGTCAGACTCGAGATTTT
CCTTTCTCTTTACCAGCTGATCCTTCTATGTATTATAAATATTAGACTCCGTAACAACTTCATTT
TTCTTTCTCTGACTTTCATGCTGCATTTTAAAAGGTATAAACAACAGAGACTGAATTAATTGAATCA
GTAATTACTTTTTCCAATACAAATTTGAATGCAGAATACTTGGAACTTGTACATGGGGAAAAATGAC
TTATCACTACATAATGTCCAATGTTTTTTCTATGTTTTGTACCAAAAATGGAAAAATATGACAATTC
CATGCAAAGAAAATGTATCTAAATATTTTTGTTAGATGATAAGAGAACTTGGTGGTCAGGACAGCAGC
CAGCTTTGTATTGTAATGCTATATTATGTAATGTGATGAAAATGTTTTGTGAAGTACTTGTAACTAA
ATTCTACAGCCATTTTCACTCCAACAGCTGTTTTTATTTTACCTCTTGGCCTAAAGTTTTTCATAA
TTTCAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_001029851.4](#)

**Summary:** The protein encoded by this gene is a cyclic nucleotide phosphodiesterase (PDE) that catalyzes the hydrolysis of the second messenger cAMP. The encoded protein, which does not hydrolyze cGMP, is resistant to several PDE inhibitors. Defects in this gene are a cause of autosomal dominant striatal degeneration (ADSD). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2010]

**Locus ID:** 8622

**MW:** 66.2