

Product datasheet for SC210618

ARHGAP30 (NM_181720) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ARHGAP30
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_181720
Insert Size:	799 bp
Insert Sequence:	>SC210618 3'UTR clone of NM_181720 The sequence shown below is from the reference sequence of NM_181720. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCTAACCTGGGAAAGGTGAAGGACTGTGATTAGGACCACAGCCCTGGGCAAAGGGACCAGCAAGTTG
TCTTGAATCTCCAGGGTTCTGACTAGCTGTCTCTCTGCAGCATGAGCAGCTGTAGTGCCCAACTCTA
TAGGCTTTGGCCCTCCAGCTTCTCTTTGACTGTGGGAGGCACTGCCTTGGTTGGTTTACCTGAACCT
GTCTCCGACACAAAGCACTTATCTCTTAGGAGATTCCCAAGAAAGTCAACAAGATCTTGTTCCAGGGA
GTGGGTCAATTGGCCAAAGGGAACATAAGGTAGGCAGAAAACCTAAAAGAGTTTGTAAAGTGAAGACTG
GAGAAATTCCTCCCTTCTCTGAGCTGTGAATCTCTCTTCATGAAAGCCAAAGGTAGAGACAGGGAGGA
CAGGGCCAGGTTAGGGCCTTCCACACACAAACACTTCTAGAGTTGCCATTCTGTTATGTTCTTGAC
CCTAAGATACCTCCTGTCCCTTTTAAATCCAGATTAAGAGAAACGTCCAGGAAGAGCTCTTGAAGCCC
TCAATATTTGTTGGAGGACTGGACTCCTCTCCAGCTCCCCACCTCTGCCTCCAGTCACCATGTGCAA
GAGAGGTCCTGTACAGATCTCTCTGGGCTCTCCTTTCTCCTTTGGAATAACTTGTTCATTTCAGGAA
AGGGAAATGGTGTCACTCAGGCCCTGGGACTGCTTCTCCAGCCAGGCTGGGGCCACAGGTCCCACTCTA
GTGAAGGTCAATGTCTCAGAATAAAAGCTGTATTTTACA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_181720.3</u>
Summary:	GTPase-activating protein (GAP) for RAC1 and RHOA, but not for CDC42.[UniProtKB/Swiss-Prot Function]
Locus ID:	257106
MW:	29.3