

Product datasheet for **SC210620**

TMEM44 (NM_138399) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Product Name: TMEM44 (NM_138399) Human 3' UTR Clone

Symbol: TMEM44

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PS100062)

ACCN: NM_138399

Insert Size: 881 bp

Insert Sequence: >SC210620 3'UTR clone of NM_138399
The sequence shown below is from the reference sequence of NM_138399. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGGACAGCACACCTGAGTGATGATGATTAACACCTTCTGGAGCCAGCTCATCAGCTCAGAGCCCAGGGT
CAGGAGTTCGTTCACTAACGCAGCGGGAATCAATCTGCACTGACACCGCGGCAGGAAGTGAAGCTGCCC
TGGCAAGTGAGGAACAGGAGCCGTCAGTGAGTGTGGCTGGGCTACATCATAGCTCATCACGGAGCTAC
GACTTTGGGTACTGCGGACAGACCTGGATAGGCCCAGCATTCTGTTCTGAAGATCACAGTTCACAGAAGC
TTTTGCTTCGTAAAGATAATCCAAAGGATCTCAGACCCGCTCTTCTTTTCCCTTCATTCCCTTGAGAG
TCAGCCATGAACGAATACCTGTAGGTTCCAGGAATGAGCTCACCTAACAGATAGCAAATGTGTCTGG
TTAGATCTCAGCAGAGCCCATTCTGCAAGACCTGGCTGAGCCAGATGAGAGGGTGGGCCCTGTGCTGGG
GGGCCTTGGGTACACACAGGAACCGAGACCTGGCTTCCACCCCCAGTCACCCACTTGGGTATCTGC
TGGAAGTATCGATAGGACTGTGTGGCCAACCAAGTGCTTGTGAGATCACTGACACTGCAAAAACAAAG
CAAAGTCTCCGGGTACAGGACTTCTCCAACCTGGCAAGGGTGTGCGCTGAGGCGGGGCTTGAGGT
GAGGGGGCTGTATGCTTCAGGAATACTAAATGCATGCAGAAGGTAAGAGGCATGATGGGAGGTGTTT
AAGCACAGCAATCCCATTTGGGAGTTATTTTGATACTGCGATGAGTAAGGGTAAGGGCGCATGGAATGG
GGCTAAGGTGGGAGTGAACACTGGGGTGAATAAATTTTAAATCAATTCAAGTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_138399.5</u>
Locus ID:	93109
MW:	32.7