

Product datasheet for **SC210621**

TMEM44 (NM_001166305) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TMEM44
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001166305
Insert Size:	881 bp
Insert Sequence:	>SC210621 3'UTR clone of NM_001166305 The sequence shown below is from the reference sequence of NM_001166305. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGGACAGCACACCTGAGTGATGATGATTAAACACCTTCTGGAGCCAGCTCATCAGCTCAGAGCCCAGGGT
CAGGAGTTCGTTCAAGTAACGCAGCGGAATCAATCTGCACTGACACCGCGGCAGGAAGTGAAGCTGCCC
TGGCAAGTGAGGAACAGGAGCCGTCACTGAGTGTGGCTGGGCTACATCATAGCTCATCACGGAGCTAC
GACTTTGGGTACTGCGGACAGACCTGGATAGGCCAGCATTCGTTCTGAAGTACAGTTTACAGAAGC
TTTTGCTTCGTAAAGATAATCCAAAGGATCTCAGACCCGCTCTTCTTTCCCTTCATTCCCTTGAGAG
TCAGCCATGAACGGAATACCTGCTAGGTTCCAGGAATGAGCTCACCTAACAGATAGCAAATGTGTCTGG
TTAGATCTCAGCAGAGCCATTCTGCAAGACCTGGCTGAGCCAGATGAGAGGGTGGGCCCTGTGCTGGG
GGGCCTTGGGTACACACAGGAACCGAGACCTGGCTTCCACCCCCAGTCACCCACTTGGGTTATCTGC
TGGAAGTTATCGATAGGACTGTGTGGCCAAACCAAGTGTGTGAGATCACTGACACTGCAAAAACAAAG
CAAACCTGCTCCGGTACCAGGACTTCTCCAACCTGGCAAGGGTGTGCGCTGAGGCGGGGCTGCAAGT
GAGGGGGCTGTATGCTTCAGGAACCTAACTAAATGCATGCAGAAGGTAAGAGGCATGATGGGAGGTGTTT
AAGCACAGCAATCCATTTGGGAGTTATTTTATACTGCGATGAGTAAGGGTAAGGGCGCATGGAATGG
GGCTAAGGTGGGAGTGAACACTGGGGTGAATAAATTTAAATCAATTCAAGTA
ACGCGTAAGCGGCCGCGCATCTAGATTTCGAAGAAAATGACCGACCAAGCGACCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001166305.2</u>
Locus ID:	93109
MW:	32.7