

Product datasheet for **SC210634**

ZNF235 (NM_004234) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF235 (NM_004234) Human 3' UTR Clone
Symbol:	ZNF235
Synonyms:	ANF270; HZF6; ZFP93; ZNF270
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004234
Insert Size:	900 bp
Insert Sequence:	>SC210634 3'UTR clone of NM_004234 The sequence shown below is from the reference sequence of NM_004234. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAGAGAGTCCACACTGGAGGGAATCTGAGAAATGAGAGGTGTGGTTCAGCCTTCAGTTAGAGCTCACA
TCTTTGTCTCCGTTGGGAAGTCCATGCTGATGATTGTGAAAAGTTCCTTTAAACTAGAAGCTTCAAAG
AATCTTGACAGGAAAGAAGTCTCTCAAATGCTGTGTTTTAGACTTGGTTAGGACATGAATCTTTACAA
ACGTCATATTTACAGATGCGAGAAAACCTTTGTTCTGTAATATGATCAGTGTTTATATCAGAGTACCA
GGTGTCCCGGTTCTCAAGATGTAACACAGCAGAAAAGTCTGTGCGAGTCTGCACAGGAGAGAGGCCTT
CCTTTAACAAAAGTATGATTGACAGAAAAAATGGATGTCAACTAAATGTAACCTCTATGTAAGGAGGT
GTTTTTCACTCTATCTCTACCATGTAGAATTGTGCTTGGCTCATAATAGGGCCAGCAATTACTGA
AGAAAGTGAAGAAAATACCTATAATTAGGACAAATACTTGAAAATTCTAGTCTGTTTTCTACTCTAAAT
TTATATTAATTTATATTCAAAGGAATCCTGCAAGTAGCCTTAATAACGTTAGTTTCAGTAAGTACTGA
AAAGTATATAATACTGCCATCCATACTAACACAAATTTGGAAAAAGACATGAATGGCTCCCATCCTTC
AGCCTTAGTTCCTTTGTACAGTGTATCATCAGTTGTGTACAGTATGATTTTTAGTGTAGTTTATAGACT
GAACAATACTGTTTTATCTAATGTGACCCTGCAATTTACATGATTAGAATTTTTGGACACTAAGCTGT
ACTGTAACCTCTGAAAGTGTAAAAGTAGTTACTGATGACAAAAAGTTTCTTAGTAAACTGAATTGTTA
ATA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_004234.4
Summary:	This gene product belongs to the zinc finger protein superfamily, members of which are regulatory proteins characterized by nucleic acid-binding zinc finger domains. The encoded protein is a member of the Kruppel family of zinc finger proteins, and contains Kruppel-associated box (KRAB) A and B domains and 15 tandemly arrayed C2H2-type zinc fingers. It is an ortholog of the mouse Zfp93 protein. This gene is located in a cluster of zinc finger genes on 19q13.2. [provided by RefSeq, Jul 2008]
Locus ID:	9310
MW:	34.3