

Product datasheet for **SC210673**

ZNF447 (ZSCAN18) (NM_001145544) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF447 (ZSCAN18) (NM_001145544) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ZSCAN18
Synonyms:	ZNF447
ACCN:	NM_001145544
Insert Size:	839 bp
Insert Sequence:	>SC210673 3' UTR clone of NM_001145544 The sequence shown below is from the reference sequence of NM_001145544. The complete sequence of this clone may contain minor differences, such as SNPs. Red =Cloning site Blue =Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

AGAGGCGCAGAGGTGAGCCGCTGTGCTGTCCCCTCCGGAGGGCCGCTTTGCCGCGGTGAATCCCAGA
CGAGGCATTGGGCCTTTCCACGCCCTGGGTGGCGGCTTCTGTGGTGTGGTGACGTCCTCTGCCTGT
GCCCTGAATCCGCTCCTGAGGCTAAGCGCTCCCAACGAGAAGGGTCCACGGGAAGCCCTCACCTGTAA
ACACACCTGGGCCAGCGCTCGCATCCGAGGGGAGCCCGGATGTGGAAGAAGACTCGGCTTTCTGCA
GCCATTTAGTGCCGCCCATGCTAGGTTATTTGACATTGTGCAGTGTAGAGTTGCCTTAAAGTGCCTGAT
CTGCCAGTGCTTTCTTCAAGTCACCTTGCCCCGATTCTCTGTTTGGCGTCCCAGGGTTGCTCAAGT
GGAAATTTGTCAGCTGTTTAGCCTTTTCGACTTGGCGTGATGTCAACTTCACTTCTAATCTGCAAAAG
CAGAAGCTGTTTCTAGTTTACCTCGCGTGTGTTTACCTATATGGAGTAGCTCGCAGAGATCACAGAAAT
GCTTGCAGCCTAAGGCAGGGTTTTAGACCGTGGGTCCCAGCCATTTAGTAAAATGGGAAATCAATTAG
CAAGTGGTCACCAGCATTACACAGCAATGAAGCAGAATAAAGTAGGCCAGAATGCATCATGTAGTAAAGG
CAAATACTGTTTTGTGAAACTTTTACCATACATCTAAATGTGAGAACTGGTTGCAATGTAAAGACATTT
CTTGCTGGGAAGTTGTGAGCAAAATAAGTTGAAAACACTAATAAAGATCTGTCTGTCTGAGCAAAGGAG

ACGCGT AAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001145544.1
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	65982