

Product datasheet for **SC210679**

ABCG1 (NM_207628) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ABCG1 (NM_207628) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ABCG1
Synonyms: ABC8; WHITE1
ACCN: NM_207628
Insert Size: 905 bp
Insert Sequence: >SC210679 3'UTR clone of NM_207628
 The sequence shown below is from the reference sequence of NM_207628. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCAGGTACAAAAATCCGGGCAGAGAGGTAAACACCTGAATGCCAGGAAACAGGAAGATTAGACACTGT
GGCCGAGGGCAGCTCTAGAATCGAGGAGGCAAGCCTGTGCCCGACCGACGACACAGAGACTTTCTGAT
CCAACCCCTAGAACCGCGTTGGGTTTGTGGGTGTCTCGTGCTCAGCCACTCTGCCAGCTGGGTTGGAT
CTTCTCTCCATTCCCCTTTCTAGCTTTAACTAGGAAGATGTAGGCAGATTGGTGGTTTTTTTTTTTTTTA
ACATACAGAAATTTAAATACCACAACTGGGGCAGAATTTAAAGCTGCAACACAGCTGGTGTGAGAGGC
TTCCTCAGTCCAGTCGCTCCTTAGCACCAGGCACCGTGGGTCTGGATGGGGAAGTGAAGCAGCCTCT
CAGCTGATGGCTGCACAGTCAGATGTCTGGTGGCAGAGAGTCCGAGCATGGAGCGATTCCATTTTATGA
CTGTTGTTTTTACATTTTCATCTTTCTAAGGTGTGTCTTTTTCCAATGAGAAGTCATTTTTGCAAGC
CAAAAGTCGATCAATCGCATTCAATTTAAGAAATTATACCTTTTTAGTACTTGCTGAAGAATGATTCAG
GGTAAATCACATACTTTGTTTAGAGAGGGCAGGGGTTAAACCGAGTCACCCAGCTGGTCTCATACATAG
ACAGCACTTGTGAAGGATTGAATGCAGGTTCCAGGTGGAGGGAAGACGTGGACACCATCTCCACTGAGC
CATGCAGACATTTTTAAAGCTATACAAAAAATTGTGAGAAGACATTGGCCAACCTTTTCAAAGTCTTT
CTTTTTCCACGTGCTTCTATTTAAGCGAAATATATTGTTTGTCTTCTCTAAAAAATAAAAAAAAAA
AAAAAAAA
ACGCGTAAGCGGCCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_207628.1
Summary:	The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. It is involved in macrophage cholesterol and phospholipids transport, and may regulate cellular lipid homeostasis in other cell types. Six alternative splice variants have been identified. [provided by RefSeq, Jul 2008]
Locus ID:	9619
MW:	34.3