

Product datasheet for **SC210696**

IRF2BPL (NM_024496) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	IRF2BPL
Synonyms:	C14orf4; EAP1; NEDAMSS
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024496
Insert Size:	889 bp

Insert Sequence: >SC210696 3'UTR clone of NM_024496
The sequence shown below is from the reference sequence of NM_024496. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GTAAAGTGAAAAAGGAGAGAGACCTTGAACCACTGGGCAGCCACCTCCTTTGCCCTAGACCAGCTCC
TCTCCAATCCTGAGGGCCCTCCCCAACCACTCGACCTCCCTCCCTCACCCCAAGGTGTAGAA
TTGTGAATATAACGAACTGCAAAAAGTTAGTCTTATGTATAGACATTATTTTCGTCGTATGTTTCTAT
ATTTTGAAACAAAGGTATGTAACCTTCTCATTTGAAGGATAAGCTGGTTGTGTTAAGCAGTATAGTAT
TGGTTGGGTCATTGTCATCATATCGTTAGCATTATTTGGTGGCAGAATGTTTGCCTAGGTACAGAATT
AATAGCCCTTAGCAACGACTGCTGCTGGTGTATTTTGTAAATGTTATGCACTCTCTGAAAGGAAAA
ACACACACAAAAGAAAAAGACTTTTTTTTTTTTTTTTTTTTTTTTGCCAAAGCCAGTGTTGCTGCCTAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAATGCTATAAAATGGTGAAAGCTTCCTTCTAAACTGC
CCCAAGTGTTGAAGTCTTCACTTTATTTTGTGTTTGTGTTTCTGTTTGTGTTTGCAGAAATGGT
AAGGGGGTGTGCGGGGGGATGGGGTGTATTTGTGCAAGTTGTGAGGGGAAATGTTTGGTTGTT
TCTACTGACCTGAATGTGTTGGATCTACACGTGTTGTTTGTGTTTGTCTTATTGATGCACGGATGCTT
TTGAACAGTAGAGCGAAATGCTAGACATGGAGAATCTGCTCTGTTGTCTTTATACATTTCTGTAGTT
AACAGAACTGTAATGTGCCTTGGAGCTTAGTAACCTGTAATAAATCAATTGATATTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024496.4</u>
Summary:	This gene encodes a transcription factor that may play a role in regulating female reproductive function. [provided by RefSeq, Jun 2012]
Locus ID:	64207
MW:	33.8