

## Product datasheet for **SC210700**

### **FXYD3 (NM\_005971) Human 3' UTR Clone**

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** FXYD3 (NM\_005971) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** FXYD3  
**Synonyms:** MAT8; PLML  
**ACCN:** NM\_005971  
**Insert Size:** 897 bp  
**Insert Sequence:** >SC210700 3'UTR clone of NM\_005971  
The sequence shown below is from the reference sequence of NM\_005971. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCATCACCCAGGCTCAGCCAAAGCTGATGAGGACAGACCAGCTGAAATTGGGTGGAGGACCGTTCT
CTGTCCCAGGTCCTGTCTCTGCACAGAACTTGAAGTCCAGGATGGAATTCTTCTCCTCTGCTGGGA
CTCCTTTGCATGGCAGGGCCTCATCTCACCTCTCGCAAGAGGGTCTCTTTGTTCAATTTTTTTAACTCT
AAAATGATTGTGCTCTGCCAAGCAGCCTGGAGACTTCCATGTGTGCATTGGGGTGGGCTTGGGGC
ACCATGAGAAGGTTGGCGTGCCTGGAGGCTGACACAGAGGCTGGCACTGAGCCTGCTTGTGGAAAA
GCCACAGGCCTGTTCCCTTGTGGCTTGGGACATGGCACAGGCCCGCCCTCTGCCTCCTCAGCCATGGG
AACCTCATATGCAATTTGGGATTTACTAGTAGCCAAAAGGAATGAAAGAGAGCTTAACCAGATGGAAC
ACTGGAACATTCCAGTGGACCCTGGACCATTCCAGGAAAAGTGGGACATAGGATCGTCCCCTATGATG
GAAGTGTTCAGACAGTTTATAATAGTAAGCCCTGTGACCCTCTCACTTACCCCGAGACCTCACTTTAT
TACAAGATCTTTCCAAATACCCAAATGTCCCTGCAAGCCCGTTAAATAATTCCCTATGCTACCCTTAAT
AACATAAATGACCACATAGTGTGAGAAGTCCAAACAAGCCTCAAAGTCCCTTGAGACTCCCAATACC
TAATAAGGCATGCGAAATGTTCTCATGAACTACCCACAACACGCCTAAAAGTCAAACACCCAAAAAT
ATCTCCTCCAATGTCCTGAAACATGAACCCAAAAAGAGACCCACAATAAACTCGTGACTTGTCCCCTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_005971.4</a>
<b>Summary:</b>	This gene belongs to a small family of FXYD-domain containing regulators of Na <sup>+</sup> /K <sup>+</sup> ATPases which share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD, and containing 7 invariant and 6 highly conserved amino acids. This gene encodes a cell membrane protein that may regulate the function of ion-pumps and ion-channels. This gene may also play a role in tumor progression. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, Oct 2008]
<b>Locus ID:</b>	5349
<b>MW:</b>	33.6