

Product datasheet for **SC210704**

BCOR (NM_001123384) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	BCOR (NM_001123384) Human 3' UTR Clone
Symbol:	BCOR
Synonyms:	ANOP2; MAA2; MCOPS2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001123384
Insert Size:	888 bp
Insert Sequence:	>SC210704 3'UTR clone of NM_001123384 The sequence shown below is from the reference sequence of NM_001123384. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGTGATCTGGCCTCAGACAACACTACTGGAGCAAGCTGGACCCACCATGTACAGTGTGTATAGTGTTA
ATCCTTGTGCATATGTGTCATAATAACAATTTCTGTAAAGAAAGGACACTATTACATATGAAAATAT
CTCTTCTTTATATAAGAGAAATTACTCCAGTCAGAAGGACTTAGAAACATGTTTTTTCTTTTAACT
TTTAAGTCAGTTTTTATGAAGTTGTTATAATGTTTCTTTACTTTTCAATGCACACATGCTTTGGGATAC
GTTTGTTTTTACTTGGAACATTTGTTTCTTTCTTTTTAAGGAGAAAAAAAATGAGTAAAGGAGCT
CCACACTTTGACTTAATTCATACAAAGCTCTGATGACAGGCCATGACTGTAGAGTGGTCAGAACTGTG
TGGTTGGTTTGAGGGAGCGAATTCGGGGAAGGCACTTGGTGATATAACTTTGTTTTGTTACAGAGTAC
CTGCTCGGGCCAGGTAATGCTATTGGATGTAATCCAGTAGTGTGTAATATAAATTCAAACCATATCCA
CACACAACAACTAATTGTATGAACTTTTATACCTAATTTAAAGCTGTGAAATTAGTTTTACGCAT
CAAACCGATTGTTTATATGTTTAAACATTTTATGCTCTTATTTAAAGAAGACTTTGAGCTATTTTTTT
CTGTACCCTGTAAATATTGAAAATAACATAATATGTTGAGGTGCTTGAAATGTACATAAACTAA
AATTTTCTGAATCGTGTGTTTATGTTTGAATCTGTGTTTAACTTTGTAAGTAAATCTCTGCCTTTG
TATTTATATTTTACAAAAATTTCTTAAAGGCAATAAACTGTTGAGGAAAGGAGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001123384.2</u>
Summary:	The protein encoded by this gene was identified as an interacting corepressor of BCL6, a POZ/zinc finger transcription repressor that is required for germinal center formation and may influence apoptosis. This protein selectively interacts with the POZ domain of BCL6, but not with eight other POZ proteins. Specific class I and II histone deacetylases (HDACs) have been shown to interact with this protein, which suggests a possible link between the two classes of HDACs. Several transcript variants encoding different isoforms have been found for this gene. A pseudogene of this gene is found on chromosome Y.[provided by RefSeq, Jun 2010]
Locus ID:	54880
MW:	35.2