

Product datasheet for **SC210724**

COMMD8 (NM_017845) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	COMMD8 (NM_017845) Human 3' UTR Clone
Symbol:	COMMD8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017845
Insert Size:	867 bp
Insert Sequence:	>SC210724 3' UTR clone of NM_017845

The sequence shown below is from the reference sequence of NM_017845. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

AGTCCTTGAAGCAGCGAATAAGGTGGTCTGCAGTTGAAA**TA**ACTGGAAATGATGAATACCAAGTCCTAT
CAGATTTTATTGCTCCAACCTTATATGGCAGAGTGAATACTGCGTGTTCCAGAAACCTTGTGATGCTTGAC
TGTTGCACCAAGGCTGAGAAAGCAGCAATATTGATATTATAAGATAAAAAATTTATCAACATTCCTTAACA
GGAAATTACATGGTTGAGAGGAAATGCATAAAATGAAAGATGAAAAATCTATAGTAGCAGTTTATATTTT
CATGATTGTTTTGCCTCATTATTAAATATTTGAGAAATCTTTGAGATACATAGTTTTATTGAAAGCTA
AAAATAGGTTCTAAAGTAATGAAAAATATAAGCACAAATATACTTGAATATTGCTTAAAGAATTGTGT
GAATAGCAACATATATTATGGATATATACTTTGTGATATTTTAAAAATAATTTTTTCAAAGAATGTAT
AAGCTGCATATATACTCAGGAGATTCCATGCTTTCTCATATTTAGAGGAAAGATTATAAAATATAAA
ATTTCTTAGAGAACACCTCTTTGTCAGAGATAAACAAGAACAATACTCTAACTTATGTGAACAGTTTT
GAGTTTATGAATTCTAGAACTAAAATCAAGAATACAGAAAAATGAAAAAACATTTTACTTCTGCGCTT
CTATGTTTGGGAAACATTGCTCTGATAAAAAATAGCTGTCATTATGCAGTGTGTATTTCAAATATGAGA
TAAGACTATGTACACATCCACTTTTGTATAAACTCAATATTGAATACTTTTGGATGTTAAATTCATTG
GAAAAACAAACCATTTGTAACCTCAGT

ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCG

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_017845.3</u>
Summary:	The protein encoded by this gene binds coiled-coil domain-containing protein 22 (CCDC22), and this complex can regulate the turnover of I-kappa-B and the activation of NF-kappa-B. [provided by RefSeq, Jul 2016]
Locus ID:	54951