

Product datasheet for **SC210728**

TCFL5 (NM_006602) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TCFL5 (NM_006602) Human 3' UTR Clone
Symbol:	TCFL5
Synonyms:	bHLHe82; CHA; E2BP-1; Figlb; SOSF1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006602
Insert Size:	1000 bp
Insert Sequence:	>SC210728 3'UTR clone of NM_006602 The sequence shown below is from the reference sequence of NM_006602. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGC
CAGAGCAGCCCCGATGGAGATCAAGTGATCGGACTGAACAGGAATCCTCGGGGGTGAACAGCCATT
CCTTCGTGACCTGTGCACGCCTTCTGCAACCCTGGAGCTCTGCTCGGCTAGTCTGACTCGAAAAGGGCG
TGACTCAAGCTGACGGGACTCCAGTAGGACTTTGAGAGCACATTTGTAAAAATATTTATCTAGACGC
AAATGCTTATCCATGAATGTCCTTTAGACCATTGGGGATGAAGCCATCTTAATAATTAGTAATAATT
AATTAGTAATAATTAGTAAGCATTTTCTCAATGCTCTGATTCCATCATGTTTTCTTAACATGATAACTT
AAAAAATTGACATCCTTTGTACTTTCTTAATCTTAAAAAGTACACGGCTTTTACTTATTTACCTTTT
AAATATGCCCCTTTAGCAATTGGAACAAGTTAAATTGTAACTAAAAACAGTTTGGAAATTTTATTTCA
TTCGTTATATCACACCCCTTGTGATGACTCTGAGTCACGTGCTGCTGTATTGCAACGTGCAGGACCAT
TTTAAACCTGTGTGCTAAAAATTTCCAGATACTTGCTTTAAAGTACTTTTGTCCACAAATGAAATAC
TGTCACAGTAGACGCTTAAATGCCACGTTTTATACCAAGAGTCATTCTTACTTCATGTGTCACAAAC
TGTGGTGTGTTGGAATTGGGTTTTCAATGAGTGGCTTTACTTATCAATCACAACAGGTAATAGCAATAG
ACGTTAGTGCAATACAAAGTCACCCTCAATAAATACTGTTAATTGGAGATGTGAGTTTGTACACAAAAC
ATCAGACTAGACCTTTGTATGGGAGAGAAATTTACTGTACATTAATTTCTTATTTTGTACTCTGTA
GCCAGTCATTAATCTGAAGGTTTAAATATATCATTTTATTGGGATGAGATCATAGCTTTACACAAATGC
TATGTAAACAAGTTACTGAATATTTTACCTCG
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006602.4</u>
Summary:	Putative transcription factor. Isoform 3 may play a role in early spermatogenesis. [UniProtKB/Swiss-Prot Function]
Locus ID:	10732
MW:	38.2