

Product datasheet for **SC210747**

TMEM156 (NM_024943) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TMEM156
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_024943
Insert Size:	966 bp

Insert Sequence: >SC210747 3'UTR clone of NM_024943
The sequence shown below is from the reference sequence of NM_024943. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAGTGCTTCCCCAATTCCAGAACTAAGTTACTTCCACAGTGCATCAGTGAGATCAATATACACGA
ATATCCCGGGCAAGTTGGACCGAGCCCTTTGAAGAATACTCAGAAGTTATTTTGTGAATGAGTAGAC
TGGAAAATGTTTGTCCAGCTGAGGATGCACAGTTGGAAAGCAGGAGGAATGCTGACTGGTTGATGAA
AACTAGCTTAAGAGCATTATTCGCTCCATGAGATCAAGGGAACAAGAGTGTTTGAAGAAGCCATTAT
GAGTCATGGAAAAAAGATGATGAAACCATGGAACAGCAAGAGAATCCCACTCTCTCTTCTTAA
AAAAAATCTATCATTATACAGCACAGAGTGGAGCCAAGTTTTTAATTTTGAAGAACCAAAAAACAGGATC
AAATATGAAAACCTTTCTTTTATTGGGCCACATTGTAGATGCTGATTTGATAATTGTTTCCTATGCAG
ATAGATTATTTTATTTACAGATTATTTAAAGGGAAGAGGCCCTGGTTGTTTATTTATATGTTTGT
TGCATTTATGAATCTTGCTGCCTTTTAGCACCAGGATGTTTTTAAAAAATTCAAAGAGGCCAGGCGCA
GTGGCTCATGCCTGTAATCCCAGCACTTTGGGATTCTGAGGTGGGAGGATCATGAGGTCAAGGGATCGA
GACCATCCTGGCCAACATGGTGAAACCTGTCTCTACTAAAAACACAAAAATTAGCTGGGTGTGGTGGT
GCGCGCCTGTAGTCCCAGCTCCTCAGGAGGCTGAGGCAGGAGAATCACTTGAACCTGGCAGGCAGAGTT
TGCAGTGAACCAAGATCACGCCACTGCATTACAGCCTGGTGCAGAGCAAGACTCTGTCTCAAAAAAAA
AAAAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAAATAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024943.3</u>
Locus ID:	80008
MW:	37.4